



**Appendix F**  
**Lavender Top-EDTA Redraw/Take Home Sample Form**

**TO BLOOD DRAWING PERSONNEL**

This blood sample is for a study sponsored by the National Institute of Health (NIH). Samples are housed at Indiana University School of Medicine. It will need to be shipped to the address below. Please use the enclosed pre-addressed FedEx Clinical Pak.

**PSDC at NCRAD**  
**Indiana University School of Medicine**  
**351 W. 10<sup>th</sup> St. TK-217**  
**Indianapolis, IN 46202**  
**Phone: 1-800-526-2839**

The kit provided contains a collection tube with which to obtain blood from the individual for research purposes. Each kit contains 1 lavender-tube and all necessary shipping supplies.

**DO NOT REFRIGERATE; STORE AT ROOM TEMPERATURE.**  
**DO NOT DRAW OR SHIP ON FRIDAY OR SATURDAY.**  
**PLEASE SHIP SAME DAY AS BLOOD IS DRAWN.**

Instructions for drawing and shipping blood samples:

1. Fill **LAVENDER TUBE** completely, if possible.
2. Invert (do not shake) tube eight to ten times after drawing blood to thoroughly mix additive with sample.
3. **Enclose this form in shipment with sample.**
4. Ship samples by **Federal Express** immediately after drawing. Use the enclosed, pre-paid Federal Express mailer. There will be no cost to you or the patient for the shipping. Consult the enclosed information sheet for packing instructions.

**KIT NUMBER (RECORDED ON LABEL):** \_\_\_\_\_

**PIDN (RECORDED ON LABEL):** \_\_\_\_\_

**STUDY SITE ID (RECORDED ON LABEL):** \_\_\_\_\_

**DATE BLOOD WAS DRAWN:** \_\_\_\_\_

**DONOR YEAR OF BIRTH:** \_\_\_\_\_

**DONOR SEX:** \_\_\_\_\_