



Longitudinal Early-onset Alzheimer's Disease Study

in collaboration with

The National Centralized Repository for Alzheimer's Disease and Related Dementias (NCRAD)

Biofluids Collection Training Slides



National Centralized Repository for
Alzheimer's Disease and Related Dementias

Contact Information

- Questions?

Please contact NCRAD Coordinators at:

- Phone: 1-800-526-2839 or 317-274-7546
- E-mail: alzstudy@iu.edu or agericks@iu.edu
- Website: www.ncrad.org

Training Overview:

- Specimen Collection Schedule
- Kit Request Module
- Specimen Labels
- Handling/Processing Study Specimens
- Sample Shipping
- NCRAD Website
- Questions?

Biofluids Collection Schedule Overview

	CI Baseline	CN Baseline	CI Month 12	CN Month 12*	CI Month 24	CN Month 24	CI Month 36	CI Month 48 /Annual visit
Serum	X	X	X	X	X	X	X	X
Plasma	X	X	X	X	X	X	X	X
DNA	X	X	X	X	X	X	X	X
Whole blood for CLIA lab testing	X							
Whole blood for long read sequencing	X	X	Collected only once over the entire course of a participant's participation in the LEADS Study. May be collected at longitudinal visits if not collected at Baseline					
RNA	X	X	X	X	X	X	X	X
PBMC	X	X	X	X	X	X	X	X
CSF	X	X	X		X	X	X	

**CN M12 CSF may be collected if CSF was not collected at Baseline.*

Biofluids Collection Schedule for CI Participants:

Sample Type	Tube Type	Number of Tubes Supplied in Kit	Aliquot Volume	Tubes to NCRAD	Ship
Whole blood for isolation of serum	Plain Red-Top Serum Blood Collection Tube (10 ml)	1	N/A	N/A	N/A
	SERUM: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	4	1.5 ml serum aliquot per 2.0 ml cryovial (red cap)	Up to 4	Frozen
Whole blood for isolation of plasma & buffy coat (for DNA extraction)	EDTA (Lavender-Top) Blood Collection Tube (10 ml)	3	N/A	N/A	N/A
	PLASMA: 2.0 ml cryovials with lavender cap (residual volume placed in 2.0 ml cryovial with blue cap)	10	1.5 ml plasma aliquot per 2.0 ml cryovial (lavender cap)	Up to 10	Frozen
	BUFFY COAT: 2.0 ml cryovial	3	1 ml buffy coat aliquot per 2.0 ml cryovial (clear cap)	3	Frozen
Whole blood for testing at the CLIA laboratory (*collected at BASELINE ONLY)	EDTA (Lavender-Top) Blood Collection tube (6ml)	1	N/A	1	Frozen
Whole blood for long read sequencing (*collected only once)	EDTA (Lavender-Top) Blood Collection tube (3ml)	1	1 ml whole blood aliquot per 2.0 ml cryovial (green cap)	3	Frozen
Whole blood for RNA extraction	PAXgene™ Blood Collection Tube (2.5 ml)	1	N/A	1	Frozen
Whole blood for PBMC	Sodium Heparin (Green-Top) Blood Collection tube (10 ml)	2	N/A	2	Ambient
CSF Collection (*not collected at 48-Month/Annual visit)	Sterile Container	Conical tubes, 15 cryovial tubes (13 orange cap, 1 blue cap, 1 yellow cap)	1.5 ml CSF aliquots per 2.0 ml cryovial (orange cap); residual volume placed in 2.0 ml cryovial with blue cap; 1-2 ml for local lab placed in 2.0 ml cryovial with yellow cap.	Up to 14	Frozen

Biofluids Collection Schedule for CN Participants:

Sample Type	Tube Type	Number of Tubes Supplied in Kit	Aliquot Volume	Tubes to NCRAD	Ship
Whole blood for isolation of serum	Plain Red-Top Serum Blood Collection Tube (10 ml)	1	N/A	N/A	N/A
	SERUM: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	4	1.5 ml serum aliquot per 2.0 ml cryovial (red cap)	Up to 4	Frozen
Whole blood for isolation of plasma & buffy coat (for DNA extraction)	EDTA (Lavender-Top) Blood Collection Tube (10 ml)	3	N/A	N/A	N/A
	PLASMA: 2.0 ml cryovials with lavender cap (residual volume placed in 2.0 ml cryovial with blue cap)	10	1.5 ml plasma aliquot per 2.0 ml cryovial (lavender cap)	Up to 10	Frozen
	BUFFY COAT: 2.0 ml cryovial	3	1 ml buffy coat aliquot per 2.0 ml cryovial (clear cap)	3	Frozen
Whole blood for long read sequencing (*collected only once)	EDTA (Lavender-Top) Blood Collection tube (3ml)	1	1 ml whole blood aliquot per 2.0 ml cryovial (green cap)	3	Frozen
Whole blood for RNA extraction	PAXgene™ Blood Collection Tube (2.5 ml)	1	N/A	1	Frozen
Whole blood for PBMC	Sodium Heparin (Green-Top) Blood Collection tube (10 ml)	2	N/A	2	Ambient
CSF Collection (*collected at BL and M24, may be collected at M12 if CSF was not collected at BL)	Sterile Containers (cryovial with yellow cap)	Conical tubes, 15 cryovial tubes (13 orange cap, 1 blue cap, 1 yellow cap)	1.5 ml CSF aliquots per 2.0 ml cryovial (orange cap); residual volume placed in 2.0 ml cryovial with blue cap; 1-2 ml for local lab placed in 2.0 ml cryovial with yellow cap.	Up to 14	Frozen

Kit Request Module

<http://kits.iu.edu/leads>



National Centralized Repository for
Alzheimer's Disease and Related Dementias

Kit Request Module

- An initial stock of kits will be delivered prior to the designated site-specific start date.
- Kits and individual supplies are available to order:
 - Blood Kit: CI Baseline (BL)
 - Blood Kit: CI Month 12 – Month 72
 - Blood Kit: CN Baseline – Month 72
 - Blood Kit: Long Read Sequencing (LRS)
 - Blood Kit: Supplemental (rarely needed)
 - CSF Kit: CI BL – M72 & CN BL – M72
 - CSF Kit: Supplemental (rarely needed)
 - LP Kit: 22G
 - LP Kit: 24G
 - Shipping Kit: Frozen Shipping
 - Shipping Kit: Ambient Blood
 - Shipping Kit: CSF (rarely needed)

Blood Kit: CI Baseline (BL)

Kit comes with the supplies necessary for the collection and processing of:

- Whole blood for isolation of serum
- Whole blood for isolation of plasma and buffy coat
- Whole blood for CLIA genetic testing
- Whole blood for isolation of RNA
- Whole blood for isolation of PBMC

Blood Kit: CI M12 – M70 & CN BL – M70

Kit comes with the supplies necessary for the collection and processing of:

- Whole blood for isolation of serum
- Whole blood for isolation of plasma and buffy coat
- Whole blood for isolation of RNA
- Whole blood for isolation of PBMC

Blood Kit: Long Read Sequencing (LRS)

Kit comes with the supplies necessary for the collection and processing of:

- Whole blood for Long Read Sequencing

CSF Kit: CI BL – M70 & CN BL – M70


Kit comes with the supplies necessary for the collection and processing of:

- Cerebral Spinal Fluid (CSF)

NCRAD Kit Request Module

1. **Reminder: only order kits/supplies that will be used within 30 days of receipt.**

1. **Enter your email address – this is how you will receive updates about your order.**



NCRAD
National Centralized Repository for
Alzheimer's Disease and Related Dementias

LEADS Kit Request System

Due to ongoing supply limitations, we ask that you please only order as many kits and extra supplies that you will be able to use in the next 30 days. Doing so allows us to fulfill as many kit requests as possible without depleting stock for other kit requests in our queue. If we are not able to fulfill any part of your request due to supplies being out of stock, we will reach out about those individually.

Please enter your email address here to receive a confirmation email after completing the survey:

* must provide value


Kit Request Module

1. Choose your site from the drop down list.
2. The coordinator name and contact information will appear.
3. Verify that this information is accurate, correct if necessary.

LEADS Site <small>* must provide value</small>	Northwestern University ▼
067: Northwestern University Cognitive Neurology and Alzheimer's Disease Center (CNADC) Northwestern University Feinberg School of Medicine % Kristine Lipowski 320 East Superior Street, Searle 12-541 Chicago, IL 60611 Phone: 312-503-2486 k-lipowski@northwestern.edu	
Is the contact name above correct? <small>* must provide value</small>	<input type="radio"/> Yes <input type="radio"/> No reset
Is the shipping address above correct? <small>* must provide value</small>	<input type="radio"/> Yes <input type="radio"/> No reset
Is the e-mail address above correct? <small>* must provide value</small>	<input type="radio"/> Yes <input type="radio"/> No reset

Study Visit Kits

CI Baseline Blood-Based Kit Qty	<input type="text" value="1"/>
CI M12 Blood-Based Kit Qty	<input type="text"/>
CI M24 Blood-Based Kit Qty	<input type="text"/>
CI M36 Blood-Based Kit Qty	<input type="text"/>
CI M48 Blood-Based Kit Qty	<input type="text"/>
CI M60 Blood-Based Kit Qty	<input type="text"/>
CN Baseline Blood-Based Kit Qty	<input type="text"/>
CN M12 Blood-Based Kit Qty	<input type="text"/>
CN M24 Blood-Based Kit Qty	<input type="text"/>
Long Read Sequencing Blood-Based Kit Qty	<input type="text"/>
LEADS Blood Supplemental Supply Kit Qty (usually need 1 at start-up)	<input type="text"/>
LEADS <u>Large</u> Frozen Blood Shipping Supply Kit Qty (will fit samples from 5 subjects)	<input type="text"/>
LEADS <u>Small</u> Frozen Blood Shipping Supply Kit Qty (will fit samples from 2 subjects)	<input type="text"/>
LEADS Ambient Shipping Supply Kit Qty (need 1 per subject)	<input type="text"/>
LEADS 22G LP Kit Qty	<input type="text"/>
LEADS 24G LP Kit Qty	<input type="text"/>
LEADS CSF Kit Qty (please provide arm and visit for each CSF kit in the "Comments" section)	<input type="text"/>
LEADS CSF Supplemental Supply Kit Qty (usually need 1 at start-up)	<input type="text"/>
LEADS CSF Shipping Supply Kit Qty (only needed when planning on shipping CSF separate from blood - rarely needed)	<input type="text"/>




- Indicate the quantity needed of each kit
- Once selected, kit components of the chosen kit will appear at the bottom of the screen (Pictured)

Each **CI (EOAD) Baseline** Blood-Based Kit Contains:

Finale IU Site: 11194

Finale All Other Sites: 11188

- 1 PAXgene™ Blood Collection Tube (2.5 ml) - CT004
- 1 Plain Red Top Serum (Red-Top) Blood Collection Tube (10 ml) - CT008
- 2 Sodium Heparin (Green-Top) Blood Collection Tube (10 ml) - CT008
- 3 EDTA (Lavender-Top) Blood Collection Tube (10 ml) - CT001
- 1 EDTA (Lavender-Top) Blood Collection Tube (6 ml) - CT003
- 9 Cryovial tube (2.0 ml) with LAVENDAR cap - CV027
- 3 Cryovial tube (2.0 ml) with RED cap - CV028
- 2 Cryovial tube (2.0 ml) with BLUE cap - CV034
- 3 Cryovial tube (2.0 ml) with CLEAR cap - CV014
- 4 Disposable graduated transfer pipette (3 ml) - CV015
- 1 50ml conical (Unwrapped) - CV0198
- 7 Bubble wrap tube sleeve - SH032 (Note to sites: Reuse 3 for frozen blood tubes for shipping)
- 1 81-cell cryobox - CV021
- 1 Resealable bag - ST002
- **Labels:**
 - 25 Pre-printed Collection and Aliquot Tube Label (1 x WBLD SERR10, 4 x SER SERR10, 3 x WBLD EDTA10, 3 x BU EDTA10, 10 x PLA EDTA10, 1 x WBLD EDTA6, 1 x WBLD RNAPXT10, & 2 x WBLD NAHEP10) - LB003
 - 5 Pre-printed Kit Number Label - LB003
 - 10 Labels for handwritten Site and LEADS ID - LB003



Do you need Extra Supplies?

* must provide value

Yes
 No



25 cell Cryobox (CV005)

5
 10

81 cell Cryobox (CV021)

5
 10

Cryovial tubes (2.0 ml) with lavender cap (CV027)

10
 25

Cryovial tubes (2.0 ml) with red cap (CV028)

10
 25

Cryovial tubes (2.0 ml) with orange cap (CV017 & CV018)

25
 50

Cryovial tubes (2.0 ml) with clear cap (CV014)

5
 10

Cryovial tubes (2.0 ml) with yellow cap (CV037)

5
 10

Cryovial tubes (2.0 ml) with blue cap (CV034)

10
 25

UPS Laboratory Paks (SH053)

5
 10

UN3373 labels (LB008)

5
 10

Biohazard label (LB009)

5

Study Visit Kits

1. Indicate if you need extra individual supplies. Selecting “yes” will cause a list of supplies to appear, and you can select which item(s) you need.

Study Visit Kits


Our standard shipping time for all orders is 3 weeks.

We can ship this kit request by: **05-10-2024**

If you need any supplies in this order prior to **05-10-2024**, you must contact the NCRAD coordinator for this study: agericks@iu.edu.

Comments

Expand



1. The module automatically calculates the date in 3 weeks. This is when you can expect your kit to ship by.
2. Click Submit.

Hints When Ordering Kits...

- For a visit where you will draw CSF, you will need an LP tray in addition to CSF kit.
- You will need 1 ambient shipping kit per blood kit.
- You will need 1 frozen shipping kit per every 4-5 subjects.
- Upon site start up, you should order 1 CSF Supplemental and Blood Supplemental kit.

NCRAD Kit Request Module: When It Must be Used

- Each site will be responsible for ordering kits (labels included) and maintaining supplies on site for scheduled participants.
- To order, sites will use the Indiana University online kit ordering module: <https://kits.iu.edu/leads>
- Allow a minimum of **3 weeks** for your order to be processed and delivered.

Specimen Labels

NCRAD



National Centralized Repository for
Alzheimer's Disease and Related Dementias

Label Type Summary

1. Kit Number Labels
2. Site and LEADS ID Labels
3. Collection and Aliquot Tube Labels
 - Differ by specimen type

Kit Number Labels



Provided by NCRAD in the kits

- Used to track patient samples and provide quality assurance
- Will be placed on the following locations:
 1. Sample and Shipment Notification Form: Blood & CSF (where applicable)
 - CSF samples will have a different kit number than the blood collection specimens
 2. Outside cryobox that houses aliquot tubes during storage and shipment
 3. Outside of the biohazard bag that houses collection tubes during ambient shipment
- Kits will include extra kit labels that should be thrown away if not used

Site and LEADS ID Label

LDS:


- Subjects will be identified by their site ID and LEADS ID
- The LEADS ID may only be available shortly before the visit
- Sites will be responsible for handwriting this onto the provided labels
 - Must use fine point permanent marker

Site and LEADS ID Label Cont.

LDS:









- Write information on label prior to adhering to tube
- Label will be placed on all collection tubes:
 - PAXgene™ Blood Collection Tube (2.5ml) for RNA
 - Plain Red Top Serum Blood Collection Tube (10ml) for Serum
 - Sodium Heparin (Green-Top) Blood Collection Tube (10ml) x 2
 - EDTA (Lavender-Top) Blood Collection Tube (10ml) for DNA and Plasma x 3
 - EDTA (Lavender-Top) Blood Collection Tube (6 ml) for CLIA lab testing ****CI Baseline ONLY****
 - EDTA (Lavender-Top) Blood Collection Tube (3ml) for LRS ****collected only once****
- Kits will include extra label

Collection and Aliquot Tube Labels



Labels to be placed on ALL collection and aliquot tubes

Collection and Aliquot Tube Labels



Every combination of Sample Type and Collection Tube that you may encounter

Look to the **Sample Type** & **Collection Tube** lines to determine what tube / cryovial the label should be placed on

Specimen Type & Collection Tube Guide

SPECIMEN TYPE ABBREVIATIONS

WBLD	-	Whole Blood
SER	-	Serum
PLA	-	Plasma
BUF	-	Buffy Coat
CSF	-	Cerebrospinal Fluid

COLLECTION TUBE ABBREVIATIONS

SERR10	10mL Serum Red-Top Tube
EDTA10	10mL EDTA Lavender-Top Tube
EDTA6	6mL EDTA Lavender-Top Tube
EDDTA3	3mL EDTA Lavender-Top Tube
RNAPXT10	10mL RNA PAXGene™ Tube
NAHEP10	10mL Sodium Heparin (NaHep) PBMC Tube
STERCNT	Sterile Container (for CSF)

Specimen Labels: Blood Collection Tubes

- All Serum, EDTA, and PAXGene™ collection tubes will have two labels:
 - Label 1: Collection Tube Label
 - Label 2: LEADS ID Label



Label 1:


	LEADS
Kit: 1000001	
1000000123	
Sample Type	
Collection Container	
	

Label 2:


LDS:
<hr/>

Please ensure the left-hand barcode is near the cap

	LEADS
Kit: 1000001	
1000000123	
WBLD	
SERR10	
	



Specimen Labels: Blood Collection Tubes



1 x 10mL Serum
Red-Top Tube

3x 10mL EDTA
Lavender-Top
Tubes

1 x 6mL EDTA
Lavender-Top
Tube

1x 3mL EDTA
Lavender-Top
Tube

1 x 10mL RNA
PAXGene
Tube


2 x 10mL Sodium
Heparin (NaHep)
PBMC Tubes

Specimen Labels: Aliquot Tubes



Cap Color	Sample Type
Red Cap	Serum
Lavender Cap	Plasma
Clear Cap	Buffy Coat
Green Cap	Whole blood
Blue Cap	Residual (plasma, serum, whole blood or CSF)
Orange Cap	CSF
Yellow Cap	CSF for local lab

- **Aliquot Tubes: 2ml cryovials**
 - Serum, Plasma, Buffy Coat, Whole Blood, and CSF
- **Collection and Aliquot tube label only**
- **Place left-hand barcode near cap**



Labeling Biologic Samples

Please...

- Label all collection and aliquot tubes before cooling, collecting, processing or freezing samples.
- Label only 1 subject's tubes at a time to avoid mix-ups.
- Wrap the label around the tube horizontally. Label position is important for all tube types.
- Make sure the label is completely adhered by rolling between your fingers.

Handling/ Processing Study Specimens



Site Required Equipment

Blood Collection/Safety Equipment

1. Personal Protective Equipment (PPE)
 - Lab Coat, Safety Glasses
2. Tourniquet
3. Alcohol Prep Pad
4. Gauze Pad
5. Butterfly Needles
6. Bandage
7. Sharps Bin and Lid

Processing/Storage Equipment

1. Centrifuge capable of ≥ 2000 rcf with refrigeration to 4°C
2. -80°C Freezer
3. Wet Ice Bucket







Blood Draw Order

Important Note

In order to ensure the highest quality samples are collected, processed, and stored, it is essential to follow the specific collection, processing, and shipment procedures detailed in the following pages. **Collection of biomarkers and CSF should be collected after a minimum 6-hour fast, preferably in the morning.** Please read the following instructions first before collecting any specimens. Have all your supplies and equipment out and prepared prior to drawing blood. **Please note that the centrifuge may take 30 minutes to cool, so please plan accordingly.** Draw blood in the following order:


1. Plain Red Top Serum Blood Collection Tube (10 ml) for Serum
2. EDTA (Lavender-Top) Blood Collection Tube (10 ml) for DNA and Plasma x 3
3. EDTA (Lavender-Top) Blood Collection Tube (6 ml) for CLIA lab testing ****CI Baseline ONLY****
4. EDTA (Lavender-Top) Blood Collection Tube (3 ml) for LRS ****only collected once per participant****
5. PAXgene™ Blood Collection Tube (2.5 ml) for RNA
6. Sodium Heparin (Green-Top) Blood Collection Tube (10 ml) x 2

Sample Collection - Blood

Tube Type	Number of Tubes Drawn	Tube Image
1. Plain Red-Top Serum Blood Collection Tube (10 ml) for Serum	x1	
2. EDTA (Lavender-Top) Blood Collection Tube (10 ml) for Plasma	x3	
3. EDTA (Lavender-Top) Blood Collection Tube (6ml) for CLIA lab testing **CI Baseline ONLY**	x1	
4. EDTA (Lavender-Top) Blood Collection Tube (3ml) for LRS **only collected once per participant**	x1	
1. PAXgene™ Blood Collection Tube (2.5 ml) for RNA	x1	
3. Sodium Heparin (Green-Top) Blood Collection Tube (10 ml) for PBMC	x2	

Aliquot Cap Colors

Cap Color	Sample Type
Red Cap	Serum
Lavender Cap	Plasma
Clear Cap	Buffy Coat
Green Cap	Whole Blood (LRS)
Blue Cap	Residual
Orange Cap	CSF
Yellow Cap	CSF for local lab



Serum Preparation (10ml Red Top Tube)

Prior to blood draw, label all tubes:


1 x 10mL Serum Tube



Label with:


LDS: _____

&




WBLD
SERR10

3 x 2ml Red Cap Cryovials & 1 x 2ml Blue Cap Cryovial



Label with:




SER
SERR10

Serum Preparation (10ml Red Top Tube)




Step One




- Store tubes at room temperature.
- Label tubes and cryovials with pre-printed subject labels prior to blood draw.

Step Two




- Collect blood in Serum Tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three



- Immediately after blood draw, invert tube 5 times to mix samples.

Step Four




- Allow blood to clot for 30 minutes.
- Within 60 minutes of blood draw, centrifuge samples at 2000 x g for 10 minutes at 4°C.

Step Five





- Must be spun, aliquoted, and stored in -80°C freezer within 2 hours of collection.



- Adhere preprinted labels to the red-cap cryovials.
- Aliquot 1.5 ml into each cryovial tube.
- If a residual aliquot is created, document specimen number and volume on Sample Notification Form.
- Store serum aliquots at -80°C until shipment.

Plain Red-Top Serum Tube (Serum Collection)



Serum Aliquots (up
to 4 possible)



Close up view
of 2.0 ml
cryovial


Plasma and Buffy Coat Preparation (10ml Lavender-Top Tube x 3)

Prior to blood draw, label all tubes:


3 x 10mL EDTA Tubes



Label with:



&




WBLD
EDTA10

9 x 2ml Purple Cap Cryovials & 1 x 2ml Blue Cap Cryovial



Label with:




PLA
EDTA10

3 x 2ml Clear Cap Cryovials




Label with:






BUF
EDTA10

Plasma and Buffy Coat Preparation (10ml Lavender-Top Tube x 3)



EDTA Tube (Plasma Collection)





Plasma
Aliquots (10
possible)



EDTA Tube (Buffy Coat Collection)

Buffy Coat
layer (mixed
with RBCs)




Buffy Coat
Aliquot
(Please use
CLEAR CAP
cryovial)

Important Note:
Buffy Coat aliquots
will be distinguished
from the plasma
aliquots through a
clear cap.


Whole Blood Preparation (6 mL Lavender-Top Tube)

Prior to blood draw, label all tubes:


1 x 6mL EDTA Tube



Label with:



&



WBLD
EDTA6

Whole Blood Preparation (6 mL Lavender-Top Tube)



Step One




- Store tubes at room temperature.
- Label tubes with pre-printed subject labels prior to blood draw.

Step Two



- Collect blood in tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three



- Immediately after blood draw, invert tube 3 times to mix sample.

Step Four



- Immediately after inversion, freeze the sample in an -80°C freezer until ready to ship.




CI Subjects at Baseline Only


Whole Blood Collection (1 x 3ml EDTA Purple Top Tube)

Prior to blood draw, label all tubes:


1 x 3mL EDTA Tube



Label with:



&




You will receive 4 labels with "WBLD" "EDTA3" and "LRS" on them. Of these 4 labels, use the label with the smallest specimen number for the 3mL EDTA tube

Labels will have an additional line to indicate they are for LRS

3 x 2ml Green Cap Cryovials & 1 x 2ml Blue Cap Cryovial



Label with:



WBLD
EDTA3

Whole Blood Collection (1 x 3ml EDTA Purple Top Tube)



Step 1




- Store tubes at room temperature.
- Labels tubes with pre-printed subject labels prior to blood draw.

Step 2




- Collect blood in EDTA Tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step 3



- Immediately after blood draw, invert tube 8-10 times to mix samples.

Step 4




- Adhere preprinted labels to the green cap cryovials.
- Aliquot 1 ml into each cryovial tube.
- If a residual aliquot is created, document specimen number and volume on Sample Notification Form.
- Store whole blood aliquots at -80°C until shipment.

Important Note: Ensure all tubes are not expired prior to collection and processing of samples.


RNA Preparation (2.5ml PAXgene™ Tube)

Prior to blood draw, label all tubes:


1 x RNA PAXGene™ Tube



Label with:



&



← WBLD
← RNAPXT10

RNA Preparation (2.5ml PAXgene™ Tube)



Step One




- Store tubes at room temperature.
- Label tubes with pre-printed labels prior to blood draw.

Step Two




- Collect blood in PAXgene™ tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three



- Immediately after blood draw, invert tubes 8-10 times to mix samples.

Step Four



- Store tubes at -80°C in a wire rack until shipment.




PBMC Preparation (10ml Sodium Heparin Tube) x 2

Prior to blood draw, label all tubes:


2 x 10mL Sodium Heparin (NaHep) PBMC Tubes



Label with:




&




WBLD

NAHEP10


PBMC Preparation (10ml Sodium Heparin Tube) x 2




- Store tubes at room temperature.
- Label tubes with pre-printed labels prior to blood draw.



- Collect blood in Sodium Heparin Tubes allowing blood to flow for 10 seconds and ensuring blood flow has stopped.



- Immediately after blood draw, invert tubes 8-10 times to mix samples.





- Store tubes at room temperature until shipment.
- Ship ambient same day as blood draw

CSF Collection and Processing

*****Important Note*****


CSF samples should be collected in the morning before breakfast and after an overnight fast. **Collection of biomarker fluids and CSF should be collected after a minimum 6-hour fast.** Only water is permitted until blood draws and the lumbar puncture are completed. Please remember to record “Last time eaten” on CSF Biological Sample and Shipment Notification Form.




CSF

Prior to blood draw, label all tubes:


1 x 50ml Sterile Container



Label with:




&




You will receive 15 labels with "CSF" and "STERCNT" on them. Of these 15 labels, use the label with the smallest specimen number for the 50ml conical.

13 x 2ml Orange Cap Cryovials & 1 x 2ml Blue Cap Cryovial



Label with:



CSF
STERCNT

*NCRAD does not provide a label for the yellow cap cryovial for CSF for your local lab.


CSF Preparation (15-20 ml total)

Step One




- Label tubes with pre-printed subject labels prior to collection.
- Pre-chill all cryovials on wet ice.

Step Two



- Collect initial 1-2ml (if bloody, collect CSF until cleared of blood) into 15 ml conical tube.
- If not bloody, transfer 1-2 ml into the yellow-cap cryovial.
- Send to local lab for testing.

Step Three




- Collect 15-20 ml total, including the 1-2 ml sent to the local lab.
- Collect sample into 2 15 ml conical tubes.

Step Four




- Place samples upright on wet ice until centrifugation begins.

Step Five




- Preferably within 15 minutes of collection, centrifuge samples at 4°C at 2000 x g for 10 minutes.

Step Six



- Using a clean transfer pipette, transfer CSF from both 15 ml conical tubes into a 50 ml conical tube, leaving the debris in the bottom.
- Gently invert the 50 ml conical tube 3-4 times to mix the sample.
- Aliquot 1.5 ml into the orange-cap cryovials.
- If a residual aliquot is created, aliquot into blue-cap cryovial. Document specimen number and volume on CSF Sample Notification Form.
- Within 2 hours of CSF collection, samples need to be spun, aliquoted and in the freezer. Store at -80°C until shipment. Record time of freezing on CSF Sample Notification Form.

Step Seven



Sample Shipping


Sample Shipment Summary

Sample Type	Processing/ Aliquoting	Tubes to NCRAD	Ship	Days to Ship
Whole blood for RNA extraction	N/A	1	Frozen	Monday-Wednesday
Whole blood (Plain Red-Top Serum Tube) for isolation of serum	1.5 ml serum aliquots per 2.0 ml cryovial (red cap) ; residual volume placed in 2.0 ml cryovial with blue cap	Up to 4	Frozen	Monday-Wednesday
Whole blood for PBMC	N/A	2	Ambient/ same day	Monday - Thursday
Whole blood (Lavender-Top EDTA) for isolation of plasma & buffy coat (for DNA extraction)	1.5 ml plasma aliquots per 2.0 ml cryovial (lavender cap) ; residual volume placed in 2.0 ml cryovial with blue cap	Up to 10	Frozen	Monday-Wednesday
	1 ml buffy coat aliquot per 2.0 ml cryovial (clear cap)	3	Frozen	Monday-Wednesday
Whole blood (Lavender-Top EDTA) for CLIA lab testing	N/A	1	Frozen	Monday-Wednesday
CSF Collection	1.5 ml CSF aliquots per 2.0 ml cryovial (orange cap); residual volume placed in 2.0 ml cryovial with blue cap; 1-2 ml for local lab placed in 2.0 ml cryovial with yellow cap.	Up to 14	Frozen	Monday-Wednesday

Ambient Sample

- Sodium Heparin/PBMC
- Only Monday-Thursday collection and same day shipping. Plan ahead to schedule UPS.
- Samples must be received at IU one day after collection.
- Do NOT draw or ship ambient samples on Friday
- Include copy of Biological Sample Shipment and Notification Form

Ambient Sample Shipping




- Place refrigerant pack in the freezer 24 hours prior to shipment.
- Place filled and labeled Sodium Heparin tubes within the slots in the absorbent pad and place in biohazard bag.
- Place the kit number label on biohazard bag.
- Place the refrigerant pack into the cooler on top of the filled biohazard bag. Place lid on cooler.
- Place the cooler in the small IATA Shipping Box.
- Place an extra copy of the “Biological Sample and Shipment Notification Form” within the shipping box along with a list of contents form.
- Close shipping box and ensure labeled with UN3373 label.
- Place box within a provided UPS ClinPak, seal, and place UPS label on outside of package.

Frozen Sample Shipping

• Ship Monday-Wednesday Only

- RNA, Serum, Plasma, Buffy Coat, CSF, Whole Blood for CLIA*, and Whole Blood for LRS* (*when applicable)
- Hold packaged samples in a -80°C freezer until pickup.
- Batch Samples together
 - 5 Cryoboxes
 - Batch shipping should be performed every 3 months or as a full shipment of specimens accumulates, whichever is sooner.

Frozen Shipping - Cryoboxes




A cryobox may have up to three kit number labels on its lid, depending on what samples were collected during the participant's visit.

This is an example of a cryobox with 3 kit number labels, one for each: regular blood kit, LRS blood kit, and CSF kit

Frozen Shipping - Cryoboxes

Aliquots from a regular Blood Kit



Please note that the yellow-cap cryovial of CSF should NOT be sent back to NCRAD. If used, you should provide it to your local lab.


Frozen Shipping - Cryoboxes



Place frozen RNA (and frozen EDTA (6ml & 3ml) tubes, when applicable) in bubble wrap tube sleeves.

Place cryobox and frozen tubes in one Biohazard Bag.

Frozen Shipping – Dry Ice Requirements



- Fully cover the cryoboxes with about 2 inches of dry ice in the provided shipper.
- Each Styrofoam shipper must contain about 45 lbs (20 kg) of dry ice.

Frozen Shipping – Dry Ice Requirements

Dry Ice label should not be covered with other stickers and must be completed or the shipping carrier will reject/return your package!



Shipping Frozen Samples

- Schedule UPS
- *Send Biological Sample and Shipment Notification Form to IU ahead of shipment*
 - *Email: alzstudy@iu.edu* or
 - *Fax: 317-321-2003*

Shipping Regulations and Training

PLEASE NOTE:

- All study personnel responsible for shipping should be certified in biospecimen shipping.
- It is the responsibility of each site to ensure that the appropriate training has been provided and conducted in regards to IATA shipping.

Please see following slides for resources.

Federal Regulations/Training

- Sites are responsible for ensuring proper training is obtained.
- Current federal and international regulations require anyone directly involved with the shipment of potentially infectious materials and other regulated biological materials (including biological specimens and cultures) **be properly trained on pertinent shipping requirements.**
 - **International Air Transport Association (IATA) Training**

<p>DGI Training Center 800-338-2291 DGItraining.com Provides IATA Certified Air Seminars and online courses</p>	<p>IATA Training Schools North America 1(514)390-6726 Europe, Africa & Middle East 41 (22) 799 2751 Asia, Australia & the Pacific 65 239 7232 www.iata.org Training schools located in 30 countries</p>
<p>Saf-T Pak Inc. www.saftpak.com Provides dangerous goods training via CD or on-site instruction for North America and Europe</p>	<p>Aiconsult Email: Airconsult@wanadoo.fr www.airconsult-bf.com</p>
<p>Bureau of Dangerous Goods LTD., TIANJIN Addr.: No.3 Yingshui road, Nankai district, Tianjin China Tel: 022-23495890 83326960 83326854 / Fax: 022-83326959 Email: cadmin@bdg-china.com.cn www.bdg-china.com.cn</p>	

UN3373 Biological Substance, Category B Training

- Biological Substance, Category B are specimens being transported for “investigational purposes”
- Recommend: investigator sites document training of category B/dangerous goods
- We recommend establishing a record of your staff’s training and date of instruction
- The training records must be made available upon request by the appropriate national authority
 - Additional information from the Department of Transportation (DOT) can be found on their website <http://hazmat.dot.gov>

Biological Sample and Shipment Notification Forms

- A copy of the sample form *must* be emailed or faxed to NCRAD prior to the date of sample arrival.
- Please include sample forms in all shipments of frozen and ambient samples.
- Email: alzstudy@iu.edu
- Fax: 317-321-2003

Biological Sample Notification Form- Blood



Participant ID: LDS _____
Biological Sample and Shipment Notification Form
 Please email or fax the form on or prior to the date of shipment



Send by E-mail or Fax prior to shipment and include a copy in each shipment.

You will need to scan and print a second hard copy so that you can include the sample form in both the ambient and frozen shipments.

To: Kelley Faber		Email: alzstudy@iu.edu		Phone: 1-800-526-2839	
General Information:				Kit #:	
From: _____				<div style="border: 1px dashed black; padding: 5px; text-align: center;">KIT BARCODE</div>	
Phone: _____					
Email: _____					
Date: _____					
Study: LEADS: <input type="checkbox"/> CI Participant <input type="checkbox"/> CN Participant				Kit # (Only if 3ml EDTA tube used):	
Visit (circle one): BASELINE M12 M24 M36 M48 M60				<div style="border: 1px dashed black; padding: 5px; text-align: center;">KIT BARCODE</div>	
Sex: <input type="checkbox"/> M <input type="checkbox"/> F					
Year of Birth: _____					
Tracking #: _____					
Blood Collection:					
1. Date Drawn (MM/DD/YYYY): _____		2. Time of Drawn (24 hour clock): _____ [HHMM]			
3. Last time subject ate (MM/DD/YYYY): _____		4. Last time subject at (24 hour clock): _____ [HHMM]			
Blood Processing:					
RNA (PAXgene Tube)			Plasma & Buffy Coat (Lavender Top Tube – 10mL)		
Total Volume of blood drawn (1 x 2.5 mL PAXgene RNA tube): _____ mL			Time spin started (24 hour clock): _____ [HHMM]		
Time PAXgene RNA tube placed in freezer (24 hour clock): _____ [HHMM]			Duration of centrifuge: _____ minutes		
Storage temperature of freezer: _____ °C			Temp of centrifuge: _____ °C		
			Rate of centrifuge: _____ xg		
			Original volume drawn (3x10 mL EDTA tubes):		
Serum (Red Top Tube)			EDTA #1: _____ mL EDTA #2: _____ mL EDTA #3: _____ mL		
Time spin started (24 hour clock): _____ [HHMM]			Time aliquoted: _____ [HHMM]		
Duration of centrifuge: _____ minutes			Number of 1.5 mL plasma aliquots created: _____		
Temp of centrifuge: _____ °C			If applicable, volume of residual plasma aliquot (less than 1.5 mL-Blue cap): _____ mL		
Rate of centrifuge: _____ x g			If applicable, specimen number of residual plasma aliquot (Last four digits): _____		
Original volume drawn (1x10 mL Serum tube): _____ mL			Time aliquots placed in freezer (24 hour clock): _____ [HHMM]		
Time aliquoted: _____ [HHMM]			Storage temperature of freezer: _____ °C		
Number of 1.5 mL serum aliquots created: _____			Buffy coat aliquot #1 (last four digits): _____		
If applicable, volume of residual serum aliquot (less than 1.5 mL-Blue cap): _____ mL			Buffy coat volume #1: _____ mL		
If applicable, specimen number of residual serum aliquot (Last four digits): _____			Buffy coat aliquot #2 (last four digits): _____		
Time aliquots placed in freezer (24 hour clock): _____ [HHMM]			Buffy coat volume #2: _____ mL		
Storage temperature of freezer: _____ °C			Buffy coat aliquot #3 (last four digits): _____		
EDTA (Lavender Top Tube – 3mL)			Buffy coat volume #3: _____ mL		
3mL EDTA tube for LRS collected? <input type="checkbox"/> Yes <input type="checkbox"/> No			Buffy coat aliquot #4 (last four digits): _____		
Original volume drawn (1x3mL EDTA tube): _____ mL			Buffy coat volume #4: _____ mL		
Time aliquoted: _____ [HHMM]			PBMC (NaHep Green Top Tube)		
Whole blood aliquot #1 (last four digits): _____			Original volume drawn (2x10mL PBMC tube): _____ mL		
Whole blood volume #1: _____ mL			EDTA (Lavender Top Tube – 6mL)		
Whole blood aliquot #2 (last four digits): _____			6mL EDTA tube for CLIA testing collected? <input type="checkbox"/> Yes <input type="checkbox"/> No		
Whole blood volume #2: _____ mL			Original volume drawn (1x6mL EDTA tube): _____ mL		
Whole blood aliquot #3 (last four digits): _____			Notes:		
Whole blood volume #3: _____ mL			_____		
Whole blood aliquot #4 (last four digits): _____			_____		
Whole blood volume #4: _____ mL			_____		

Biological Sample Notification Form- Blood



Participant ID: LDS _____
Biological Sample and Shipment Notification Form
 Please email or fax the form on or prior to the date of shipment



There should be a second Kit Number Label placed on the form if a 3mL EDTA tube was collected and aliquoted for LRS

To: Kelley Faber		Email: alzstudy@iu.edu		Phone: 1-800-526-2839	
General Information:				Kit #:	
From: _____				<div style="border: 1px dashed black; padding: 5px; text-align: center;">KIT BARCODE</div>	
Phone: _____					
Email: _____					
Date: _____					
Study: LEADS: <input type="checkbox"/> CI Participant <input type="checkbox"/> CN Participant				<div style="border: 1px dashed black; padding: 5px; text-align: center;">KIT BARCODE</div>	
Visit (circle one): BASELINE M12 M24 M36 M48 M60					
Sex: <input type="checkbox"/> M <input type="checkbox"/> F					
Year of Birth: _____					
Tracking #: _____					
Blood Collection:					
1. Date Drawn (MM/DD/YYYY): _____		2. Time of Drawn (24 hour clock): _____ [HHMM]			
3. Last time subject ate (MM/DD/YYYY): _____		4. Last time subject at (24 hour clock): _____ [HHMM]			
Blood Processing:					
RNA (PAXgene Tube)			Plasma & Buffy Coat (Lavender Top Tube – 10mL)		
Total Volume of blood drawn (1 x 2.5 mL PAXgene RNA tube): _____ mL			Time spin started (24 hour clock): _____ [HHMM]		
Time PAXgene RNA tube placed in freezer (24 hour clock): _____ [HHMM]			Duration of centrifuge: _____ minutes		
Storage temperature of freezer: _____ °C			Temp of centrifuge: _____ °C		
			Rate of centrifuge: _____ xg		
			Original volume drawn (3x10 mL EDTA tubes): _____ mL		
Serum (Red Top Tube)			EDTA #1: _____ mL EDTA #2: _____ mL EDTA #3: _____ mL		
Time spin started (24 hour clock): _____ [HHMM]			Time aliquoted: _____ [HHMM]		
Duration of centrifuge: _____ minutes			Number of 1.5 mL plasma aliquots created: _____		
Temp of centrifuge: _____ °C			If applicable, volume of residual plasma aliquot (less than 1.5 mL-Blue cap): _____ mL		
Rate of centrifuge: _____ xg			If applicable, specimen number of residual plasma aliquot (Last four digits): _____		
Original volume drawn (1x10 mL Serum tube): _____ mL			Time aliquots placed in freezer (24 hour clock): _____ [HHMM]		
Time aliquoted: _____ [HHMM]			Storage temperature of freezer: _____ °C		
Number of 1.5 mL serum aliquots created: _____			Buffy coat aliquot #1 (last four digits): _____		
If applicable, volume of residual serum aliquot (less than 1.5 mL-Blue cap): _____ mL			Buffy coat volume #1: _____ mL		
If applicable, specimen number of residual serum aliquot (Last four digits): _____			Buffy coat aliquot #2 (last four digits): _____		
Time aliquots placed in freezer (24 hour clock): _____ [HHMM]			Buffy coat volume #2: _____ mL		
Storage temperature of freezer: _____ °C			Buffy coat aliquot #3 (last four digits): _____		
EDTA (Lavender Top Tube – 3mL)			Buffy coat volume #3: _____ mL		
3mL EDTA tube for LRS collected? <input type="checkbox"/> Yes <input type="checkbox"/> No			Buffy coat aliquot #4 (last four digits): _____		
Original volume drawn (1x3mL EDTA tube): _____ mL			Buffy coat volume #4: _____ mL		
Time aliquoted: _____ [HHMM]			PBMC (NaHep Green Top Tube)		
Whole blood aliquot #1 (last four digits): _____			Original volume drawn (2x10mL PBMC tube): _____ mL		
Whole blood volume #1: _____ mL			EDTA (Lavender Top Tube – 6mL)		
Whole blood aliquot #2 (last four digits): _____			6mL EDTA tube for CLIA testing collected? <input type="checkbox"/> Yes <input type="checkbox"/> No		
Whole blood volume #2: _____ mL			Original volume drawn (1x6mL EDTA tube): _____ mL		
Whole blood aliquot #3 (last four digits): _____			Notes:		
Whole blood volume #3: _____ mL			_____		
Whole blood aliquot #4 (last four digits): _____			_____		
Whole blood volume #4: _____ mL			_____		

Biological Sample Notification Form – CSF

Send by E-mail or Fax prior to shipment and include a copy in each shipment



Participant ID: LDS



CSF Sample and Shipment Notification Form

Please email or fax the form on or prior to the date of shipment.

To: Kelley Faber		Email: alzstudy@iu.edu		FAX: 317-321-2003		Phone: 1-800-526-2839	
General Information:							
From: <input type="text"/>		Date: <input type="text"/>		Phone: <input type="text"/>		Email: <input type="text"/>	
Study: LEADS <input type="checkbox"/> CI Participant		<input type="checkbox"/> CN Participant		Kit #:		KIT BARCODE	
Visit (circle one): BASELINE M12 M24 M36 M48							
Sex: <input type="checkbox"/> M <input type="checkbox"/> F		Year of Birth: <input type="text"/>		CSF Collected? <input type="checkbox"/> Yes <input type="checkbox"/> No			
Tracking #: <input type="text"/>		Gauge needle used for LP: <input type="checkbox"/> 22G <input type="checkbox"/> 24G					
CSF Collection:							
1. Date of Collection: <input type="text"/>		2. Time of Collection: 24 hour clock: <input type="text"/> [HHMM]		3. Last time subject ate: Date: <input type="text"/>		4. Last time subject ate: 24 hour clock: <input type="text"/> [HHMM]	
5. Collection process: <input type="checkbox"/> Gravity Method OR <input type="checkbox"/> Aspiration							
CSF Processing:							
Time spin started: 24 hour clock: <input type="text"/>		<input type="text"/> [HHMM]		Duration of centrifuge: <input type="text"/>		minutes	
Temp of centrifuge: <input type="text"/> °C		Rate of centrifuge: <input type="text"/>				x g	
Total amount of CSF collected (mL): <input type="text"/>		<input type="text"/>				mL	
Time aliquoted: <input type="text"/>		<input type="text"/>				[HHMM]	
Number of 1.5 mL aliquots created (up to 14 total): (Orange cap cryovials): <input type="text"/>		<input type="text"/>				x 1.5 mL	
If applicable, volume of CSF residual aliquot (less than 1.5 mL): (Blue cap cryovial): <input type="text"/>		<input type="text"/>				mL	
If applicable, specimen number of residual aliquot tube: (Last four digits): <input type="text"/>		<input type="text"/>					
Time frozen: <input type="text"/>		<input type="text"/>				[HHMM]	
Storage temperature of freezer: <input type="text"/>		<input type="text"/>				°C	
Notes: <input type="text"/>							

NCRAD Website

Helpful Pages

- <https://ncrad.iu.edu/contact/holiday-closures>
<https://ncrad.iu.edu/contact/friday-blood-draws>

[Home](#) / [Contact](#) / [Friday Blood Draws](#)

WHAT TO DO FOR FRIDAY BLOOD DRAWS

NCRAD is not open for business on Saturday or Sunday; therefore, we ask that no samples be shipped on a Friday. We cannot guarantee the conditions in which the samples will be held by the shipping courier over the weekend. It is important to have plans in place for each type of sample to be held over the weekend prior to shipping. Please refer to the table below for how to handle samples drawn on a Friday.

When possible, please only ship frozen samples on Monday-Wednesday. There is always the potential for an unexpected shipping courier delay and by shipping Monday through Wednesday there should be enough time to receive the samples before the weekend.

SAMPLE TYPE	TUBE TYPE	PRODUCT	SHIPMENT METHOD	FRIDAY DRAW INSTRUCTIONS
Whole Blood	Sodium Heparin	PBMC	Ambient	DO NOT DRAW ON FRIDAY. Must be drawn on Monday – Thursday.
Whole Blood	EDTA Tube	DNA Only	Ambient	Do NOT refrigerate. Please keep sample at room temperature until specimen can be shipped via next day delivery methods the following Monday.
Whole	EDTA Tube	DNA Only	Frozen	Whole blood in EDTA may be frozen in a -80°C freezer within 5 days of collection and shipped frozen on dry ice to NCRAD to remain within

[Home](#) / [Contact](#) / [Holiday Closures](#)

HOLIDAY CLOSURES

DATE	HOLIDAY
January 1	New Year's Day
3 rd Monday in January	Martin Luther King, Jr Day
4 th Monday in May	Memorial Day
June 19	Juneteenth (observed)
July 4	Independence Day (observed)
1 st Monday in September	Labor Day
4 th Thursday in November	Thanksgiving
4 th Friday in November	Friday after Thanksgiving
December 25	Christmas

LEADS Active Study Page



Donate

Get Updates

- ABOUT
- BANK SAMPLES
- ACCESS SAMPLES
- BIOMARKER ANALYSIS
- COORDINATE STUDIES
- GENETICS & FAMILIES
- CONTACT

COORDINATE STUDIES

LEADS

- NIA-AD FBS
- NAPS2
- PATH
- APOE
- 4RTNI-2
- 90+ Study
- ABC-DS
- ACAD
- ACE
- ACT
- ADCFB
- ADNI-3,4
- AGMP
- ALLFTD
- BBBSR
- BenfoTeam
- BEYONDD
- CADASIL
- CLARITI

Home / Coordinate Studies / LEADS

LEADS ACTIVE STUDY PAGE

Welcome LEADS Study staff, coordinators, and PI's.

This section encompasses study specific tools and videos for your reference. If you have any questions, comments, or new ideas please contact NCRAD by **email** or phone **1-800-526-2839** or directly at **317-278-8413**.



CI (COGNITIVELY IMPAIRED) PARTICIPANTS					
	BASELINE	M12	M24	M36	M48
RNA	✓	✓	✓	✓	✓
Serum	✓	✓	✓	✓	✓
Plasma	✓	✓	✓	✓	✓
DNA	✓	✓	✓	✓	✓
PBMC	✓	✓	✓	✓	✓
EDTA for CLIA testing**	✓				
CSF	✓	✓	✓	✓	

**Please note that this EDTA tube is used for the purpose of confirmation testing and is only drawn during the Baseline visit for CI participants.



Contact Information

- Questions?

Please contact NCRAD Coordinators at:

- Phone: 1-800-526-2839 or 317-278-1133
- E-mail: alzstudy@iu.edu or agericks@iu.edu
- Website: www.ncrad.org