

ARTFL-LEFFTDS Manual Update Training: V 5.2016



NCRAD

Training Overview:

- Study Update Overview
- Kit/Kit Request Module Update
- Specimen Label Update
- Sample Handling/Processing Update
- Sample Shipping Update
- Sample Form Updates
- Questions?

Study Update Overview – Top 3 Changes





1. Addition of new biospecimens to both ARTFL and LEFFTDS participants
 - Additional EDTA Tubes and Serum Tube
2. New blood draw order
 - PAXgene™ Tube now drawn LAST
3. Shipping supplies have been updated for the increase in specimens collected
 - Introduction of 81-slot cryobox

Updated ARTFL and LEFFTDS Study Specimens

	ARTFL	LEFFTDS (Visit 1)	LEFFTDS (Visit 2)	LEFFTDS (Visit 3)
DNA (Buffy Coat)	X	X	X	X
Plasma	X	X	X	X
PBMC	X	X	X	X
Serum	X	X	X	X
RNA	X	X	X	X
CSF	X*	X	X	X

*Select participant population Only

New Blood Draw Order

Tube Type	Number of Tubes Drawn	Tube Image
1. EDTA (Lavender-Top) Tube (10 ml)	x3	 A clear plastic 10 ml EDTA tube with a purple cap. The label is partially visible, showing 'P366543', 'BA (V2)', 'SERIAL', '47149', and '21507'.
2. Sodium Heparin (Green-Top) Tube (10 ml)	x1	 A clear plastic 10 ml Sodium Heparin tube with a green cap. The label is white with green text and includes 'BD Vacutainer Sodium Heparin', '10 ml Unit', '473114', '4134502', and '2016-05 EXPIRES'.
3. Serum Determination (Red-Top) Tube (10 ml)	x1	 A clear plastic 10 ml Serum Determination tube with a red cap. The label is white with red text and includes 'BD Vacutainer Serum Separator', '10 ml Unit', '473114', '4134502', and '2016-05 EXPIRES'.
4. PAXgene™ Tube (2.5 ml)	x3	 A clear plastic 2.5 ml PAXgene tube with a red cap. The label is white with red and black text and includes 'PAXgene', '2.5 ml', 'PAXgene', '5231724', and '2017-05-31'.

Kit Request Module Update

<http://kits.iu.edu/artfl-lefftds/>



NCRAD

Available Kits:


1. ARTFL-LEFFTDS Universal Blood Kit (V 10.2015)
2. ARTFL-LEFFTDS Universal Blood Kit (V 5.2016)
3. Optional NCRAD Kit - CSF
4. Green Top-Sodium Heparin Tube Redraw/Take Home Kit
5. Lavender Top-EDTA Tube Redraw/Take Home Kit
6. Frozen Shipping Supply Kit (V 10.2015)
7. Frozen Shipping Supply Kit (V 5.2016)
8. ARTFL/LEFFTDS Supplemental Kit
9. Individual Supplies

Please note: ARTFL Base and LEFFTDS Cycle 1, Visit 1 kits will no longer be available.

ARTFL-LEFFTDS Universal Blood Kit

Quantity	ARTFL-LEFFTDS Universal Blood Kit Components
3	EDTA (Lavender-Top) Blood Collection Tube (10 ml)
1	Sodium Heparin (Green-Top) Blood Collection Tube (10 ml)
1	Serum Determination (Red-Top) Blood Collection Tube (10 ml)
3	PAXgene™ Blood Collection Tube (2.5 ml)
30	Cryovial tube (2 ml) with lavender cap
10	Cryovial tube (2 ml) with red cap
3	Cryovial tube (2 ml) with clear cap
2	Cryovial tube (2 ml) with blue cap
4	Disposable graduated transfer pipette
51	Pre-printed labels for blood collection and aliquot tubes
5	Pre-printed labels with kit number
9	Labels for handwritten Site and RAVE ID
1	Cryovial tube box (holds up to 81 cryovials)
1	Shipping Supplies for ambient shipment of PBMC: Small canister (95kPa Certified) Small IATA shipping box Absorbent tube sleeve Sheet of bubble wrap List of contents Fed Ex Clinic Pak FedEx return airbill


Automatically includes all shipping supplies to safely ship PBMC sample to NCRAD




Specimen Label Update






Three Label Types



Kit Number



Site and RAVE ID



Collection and Aliquot Tube

Collection Tubes - BLOOD



EDTA
Tube



Sodium Heparin
Tube



Serum Determination
Tube







PAXgene™
Tube

Handling/Processing Study Specimens

- Updated Plasma and Buffy Coat Processing Schematic
- NEW Serum Processing Schematic



New Blood Draw Order

Tube Type	Number of Tubes Drawn	Tube Image
1. EDTA (Lavender-Top) Tube (10 ml)	x3	
2. Sodium Heparin (Green-Top) Tube (10 ml)	x1	
3. Serum Determination (Red-Top) Tube (10 ml)	x1	
4. PAXgene™ Tube (2.5 ml)	x3	

Plasma and Buffy Coat Preparation (10 ml Lavender-Top Tube)



Step One



- Store tubes at room temperature.
- Label tubes with pre-printed subject labels prior to blood draw.

Step Two




- Collect blood in EDTA Tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three




- Immediately after blood draw, invert tubes 8-10 times to mix samples.

Step Four



- Centrifuge samples at 1500 x g for 15 minutes at 4°C.
- EDTA Tubes need to be spun, aliquoted, and in the freezer within 2 hours from the time of collection.

Step Five



- Label cryovial tubes with preprinted labels.
- Aliquot 0.5 ml into each cryovial tube.
- If residual aliquot is created, use blue cap to indicate volume difference and document Specimen Number on Biological Sample and Shipment Notification Form.
- Store plasma aliquots at -80°C until shipment.

Step Six



- Label cryovial tube with preprinted label.
- Using a clean transfer pipette, collect the buffy coat (may have residual plasma and some RBCs included).
- Transfer the buffy coat into the cryovial tube.
- Store buffy coat aliquot at -80°C until shipment.

- ❖ **Three EDTA (Lavender-Top) Blood Collection Tubes will yield approximately:**
 - ❖ **24 to 30 Plasma aliquots with 0.5 ml per aliquot**
 - ❖ **3 Buffy Coat aliquots with 1.0 ml per aliquot**



New Blood Collection Tube: Serum Determination Tube



Serum Determination
Tube (unfilled)

Serum Determination Tube
(Immediately after blood draw – pictured below)

** Please note: After standing at room temperature for
30 minutes, blood will be clotted and immobile**




Serum Determination
Tube (after centrifuged)

Serum Preparation (10 ml Red-Top Tube)




Step One




- Store tubes at room temperature.
- Label tubes with pre-printed subject labels prior to blood draw.

Step Two



- Collect blood in Serum Tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three




- Immediately after blood draw, invert tubes 8-10 times to mix samples.

Step Four




- Allow blood to clot for 30 minutes.
- Centrifuge samples at 1500 x g for 15 minutes at 4°C.
- Serum samples need to be spun, aliquoted, and in the freezer within 2 hours from the time of collection.

Step Five



- Label cryovial tubes with preprinted labels.
- Aliquot 0.5 ml into each cryovial tube.
- If residual aliquot is created, use blue cap to indicate volume difference and document Specimen Number on Biological Sample and Shipment Notification Form.
- Store serum aliquots at -80°C until shipment.









Sample Shipments/ Shipping Supplies




Additional Kit Supplies

- Cryovial Cap Colors:

Cap Color	Sample Type	Cap Image
Lavender	Plasma	
Clear	Buffy Coat	
Red	Serum	
Clear	CSF Aliquot (0.5 ml)	
Orange	CSF Aliquot (1.0 ml) and CSF Aliquot to local lab	
Blue	Residual Aliquot (Plasma, Serum, or CSF)	

Additional Kit Supplies Cont.

- 81-Slot Cryobox and Large Biohazard Bag



Additional Kit Supplies Cont.

- Large biohazard bag to accommodate 81-Slot Cryobox
- Frozen shipments should now only include biospecimens for 4 participants



Sample Shipment Summary

Sample Type	ARTFL	LEFFTDS	Processing/ Aliquoting	Tubes to NCRAD	Ship
Whole blood (Lavender-Top EDTA) for isolation of plasma & buffy coat (for DNA extraction)	Yes (V1)	Yes (V1, V2, V3)	0.5 ml plasma aliquots per 2 ml cryovials	24-31	Frozen
	Yes (V1)	Yes (V1, V2, V3)	1 ml buffy coat aliquot per 2 ml cryovial	3	Frozen
Whole blood (Green-Top Sodium Heparin) for isolation of PBMCs	Yes (V1)	Yes (V1, V2, V3)	N/A	1	Ambient
Whole blood (Red- Top Serum) for isolation of serum	Yes (V1)	Yes (V1, V2, V3)	0.5 ml serum aliquots per 2 ml cryovials	8-11	Frozen
Whole blood (PAXgene™) for RNA extraction	Yes (V1)	Yes (V1, V2, V3)	N/A	3	Frozen
CSF	Some	Yes (V1, V2, V3)	0.5 ml and 1 ml CSF aliquots per 2 ml cryovials	Up to 33	Frozen

Ambient Sample Shipment

- Ship the sample(s) to IU on the day of collection.
- ***Sample(s) must be received at IU a day after collection.***
- ***Only ship Monday through Thursday***
 - ***No Saturday Shipments***

Frozen Sample Shipment

- All other samples
 - Plasma, Buffy Coat, Serum, PAXgene™, and CSF
 - **Ship Monday-Wednesday Only**
- Hold packaged samples in a -80°C freezer until pickup.
- Batch Samples together
 - 4 participant samples (Plasma, Buffy, Serum, CSF, RNA)
 - **Batch shipping should be performed quarterly or as a full shipment of specimens accumulates, whichever is sooner.**

Biological Sample and Shipment Notification Form Update

Appendix B

Biological Sample and Shipment Notification Form

Please email or fax the form on or prior to the date of shipment.

To: Kelley Faber Email: alzstudy@iu.edu FAX: 317-278-1100 Phone: 1-800-526-2839

General Information: FedEx tracking #: _____

From: _____ Site: _____

Phone: _____ Fax: _____

Email: _____ Date: _____

Study: LEFFTDS ARTFL LEFFTDS and ARTFL Kit #: _____

Visit: _____

Site ID: _____ RAVE #: _____

Fam #: _____ Sex: M F Year of Birth: _____ CSF Sample Donated? Yes No

Blood Collection:

1. Date Drawn: _____	2. Time of Draw: _____ <input type="checkbox"/> AM <input type="checkbox"/> PM
3. Last time subject ate: Date: _____	4. Last time subject ate: Time: _____ <input type="checkbox"/> AM <input type="checkbox"/> PM
5. Sodium heparin tube (PBMC) drawn: <input type="checkbox"/> Yes <input type="checkbox"/> No	
6. Total volume of blood drawn into 3 x 2.5 ml PAXgene RNA tubes: _____ mL	
• Were the PAXgene™ tubes the last tubes drawn? <input type="checkbox"/> Yes <input type="checkbox"/> No	

Blood Processing:

Plasma (EDTA Tube)		Serum (Serum Determination Tube)	
Time spin started:	_____ <input type="checkbox"/> AM <input type="checkbox"/> PM	Time spin started (within 30 minutes of draw time):	_____ <input type="checkbox"/> AM <input type="checkbox"/> PM
Original volume drawn (3x10 mL EDTA tube):	_____ mL	Original volume drawn (1x10 mL Serum tube):	_____ mL
Number of 0.5 mL plasma aliquots created (24-30 total): (Lavender cap cryovial):	_____ x 0.5 mL	Number of 0.5 mL serum aliquots created (8-10 total): (Red cap cryovial):	_____ x 0.5 mL
If applicable, volume of residual plasma aliquot (less than 0.5 mL): (Blue cap cryovial):	_____ mL	If applicable, volume of residual serum aliquot (less than 0.5 mL): (Blue cap cryovial):	_____ mL
If applicable, specimen number of residual aliquot: (Last four digits)	_____	If applicable, specimen number of residual aliquot: (Last four digits)	_____
Buffy coat aliquots created (one per EDTA tube): (Clear cap cryovial):	_____		
Time aliquots placed in freezer:	_____ <input type="checkbox"/> AM <input type="checkbox"/> PM	Time aliquots placed in freezer:	_____ <input type="checkbox"/> AM <input type="checkbox"/> PM

Notes: _____

- ❖ Now includes expanded blood processing section for both **plasma** and **serum**, new field for Visit Number, and expanded date/time of last time subject ate.
- ❖ Both ARTFL and LEFFTDS participants should have all aspects of form completed prior to submission

CSF Sample and Shipment Notification Form Update

Appendix C

CSF Sample and Shipment Notification Form

Please email or fax the form on or prior to the date of shipment.

To: Kelley Faber Email: alzstudy@iu.edu FAX: 317-278-1100 Phone: 1-800-526-2839	
General Information:	
FedEx tracking #: _____	
From: _____	Site: _____
Phone: _____	Fax: _____
Email: _____	Date: _____
Study: <input type="checkbox"/> LEFFTDS <input type="checkbox"/> ARTFL <input type="checkbox"/> LEFFTDS and ARTFL	Kit #: <div style="border: 1px dashed black; padding: 5px; text-align: center;">KIT BARCODE</div>
Site ID: _____ RAVE #: _____	Visit: _____
Fam #: _____ Sex: <input type="checkbox"/> M <input type="checkbox"/> F Year of Birth: _____	
CSF Collection:	
1. Date of Collection: _____	2. Time of collection: _____ <input type="checkbox"/> AM <input type="checkbox"/> PM
3. Last time subject ate: Date: _____	4. Last time subject ate: Time: _____ <input type="checkbox"/> AM <input type="checkbox"/> PM
5. Collection process: <input type="checkbox"/> Gravitational OR <input type="checkbox"/> Pull	
CSF Processing:	
1. Total number of CSF aliquot tubes: _____	4. Total number of CSF 0.5 mL aliquots transferred to NCRAD: (Clear cap cryovial): _____
2. Total amount of CSF collected (mL): _____	5. Total number of CSF 1.0 mL aliquots transferred to NCRAD: (Orange cap cryovial): _____
3. Time frozen: _____ <input type="checkbox"/> AM <input type="checkbox"/> PM	6. If applicable, volume of CSF residual aliquot (less than 0.5 mL): (Blue cap cryovial): _____
7. If applicable, specimen number of residual aliquot tube: (Last four digits): _____	
Notes: _____	

- ❖ Now includes Visit Number
- ❖ Please note that fasting information is obtained through “Last Time Subject Ate”
- ❖ Both ARTFL and LEFFTDS participants should have all aspects of form completed prior to submission

Questions?

- Please utilize the “Chat” window at the bottom of the Adobe Connect Screen.

Or

- Contact Ashley Vetor at:
 - Phone: 317-278-9546
 - Email: abozell@iu.edu