



VIVAMIND MOP Protocol Update: V8.23.2023

Section	Summary of Changes
3.3	Added Juneteenth to the list of holidays observed by Indiana University
7.1	Added Sarstedt CSF Tube label
7.1	Added picture of 2ml Sarstedt Tube pictured with the Sarstedt CSF Tubel label applied
8.2	Updated General CSF Collection Methods and CSF Collection schematic
10.1	Updated Frozen Shipping instructions and clarified that pelleted dry ice should used for all NCRAD frozen shipments
12.0	Updated CSF Appendix C highlighting that the 2ml Sarstedt MUST be collected

VIVA MIND

Manual of Procedures

National Centralized Repository for Alzheimer's Disease and Related Dementias (NCRAD)

A Trial to Evaluate the Efficacy and Safety of PQ912 in Patients with Early AD (VIVA-MIND) –

Phase 2A

Biospecimen Collection, Processing, and Shipment Manual

Version 8.23.2023



Table of Contents

1.0	Abbreviations	4
2.0	Purpose.....	4
3.0	NCRAD Information	5
3.1	NCRAD Contacts	5
3.2	Hours of Operation	5
3.3	Holiday Observations	6
4.0	Globally Unique Identifier (GUID).....	6
5.0	NCRAD Laboratory Information.....	7
5.1	Site Required Equipment	7
5.2	Biospecimens Sent to NCRAD	8
5.3	Biospecimen Collection Charts.....	10
5.3.1	Biospecimen Collection for Screening Visit.....	10
5.3.2	Biospecimen Collection for Baseline Pre-dose, Weeks 4/8/16/48 Pre-dose and Post-dose Visits.....	11
5.3.3	Biospecimen Collection for Weeks 24 and 72 Pre-dose and Post-dose Visit.....	12
5.3.4	Biospecimen Collection for Early Termination Pre-dose and Post-dose Visit.....	13
6.0	Specimen Collection Kits, Shipping Kits, and Supplies.....	14
6.1	Specimen Collection Kit Contents.....	14
6.2	Kit Supply to Study Sites	17
7.0	Blood Collection and Processing Procedures	18
7.1	Labeling Samples	18
7.2	Video List	20
7.3	Filling Aliquot Tubes (Plasma, Serum, and CSF).....	20
7.4	Serum (Red-top) Tube (10 mL) for Serum	22
7.5	EDTA (Lavender-Top) Blood Collection Tube (10 mL) for Plasma and Buffy Coat	25
8.0	Cerebrospinal Fluid Collection and Processing Procedures.....	30
8.1	Scheduling the LP.....	30
8.2	Performing the LP	30
9.0	Sample Redraws	35
10.0	Packaging and Shipping Instructions	36
10.1	Frozen Shipping Instructions.....	37
11.0	Data Queries and Sample Reconciliation	43
12.0	Appendices	44

1.0 Abbreviations

AD	Alzheimer's Disease
ADCS	Alzheimer's Disease Cooperative Study
BL	Baseline Visit
CSF	Cerebrospinal Fluid
DNA	Deoxyribonucleic Acid
EDTA	Ethylene Diamine Tetra-acetic Acid
ET	Early Termination Visit
GUID	Globally Unique Identifier
IATA	International Air Transport Association
IUGB	Indiana University Genetics Biobank
LP	Lumbar Puncture
NCRAD	National Centralized Repository for Alzheimer's Disease and Related Dementias
PHI	Protective Health Information
PK	Pharmacokinetics
RBCs	Red Blood Cells
RCF	Relative Centrifugal Force
RPM	Revolutions Per Minute

2.0 Purpose

The purpose of this manual is to provide VIVA-MIND staff (PIs, study coordinators, and the sample collection and processing teams) at the various study sites with instructions for collection and submission of biological samples for VIVA-MIND study visits. It includes instructions for biospecimen submission to the National Centralized Repository for Alzheimer's Disease and Related Dementias (NCRAD) located at Indiana University. The following samples may be collected at each study visit:

- Serum
- Plasma
- Buffy Coat (DNA Extraction)
- CSF

This manual includes instructions for collection of blood and CSF, fractionation of blood from collection tubes, aliquoting, labeling, storage prior to shipping, and shipping to NCRAD.

These procedures are relevant to all study personnel responsible for processing blood specimens to be submitted to NCRAD for the VIVA-MIND protocols.



3.0 NCRAD Information

3.1 NCRAD Contacts

Tatiana Foroud, PhD, Core Leader

Email: tforoud@iu.edu

Kelley Faber, MS, CCRC, Project Manager

Phone: 317-274-7360

Email: kelfaber@iu.edu

Diont'e Keys, B.S., CCRC, Clinical Research Coordinator

Phone: 317-274-7890

Email: dlkeys@iupui.edu

General NCRAD Contact Information

Phone: 1-800-526-2839

Email: alzstudy@iu.edu

Website: www.ncrad.org

VIVA-MIND Study Specific Webpage: <https://ncrad.org/resource/viva-mind.html>

Sample Shipment Mailing Address

NCRAD

Indiana University School of Medicine

351 W. 10th St

TK-217

Indianapolis, IN 46202

3.2 Hours of Operation

Indiana University business hours are from 8 AM to 5 PM Eastern Time, Monday through Friday.

Frozen samples must be shipped **Monday-Wednesday only**.

Check weather report to make sure impending weather events (blizzards, hurricanes, etc.) will not affect the shipping or delivery of the samples.

3.3 Holiday Observations

Date	Holiday
January 1	New Year's Day
3 rd Monday in January	Martin Luther King, Jr Day
4 th Monday in May	Memorial Day
June 19	Juneteenth (observed)
July 4	Independence Day (observed)
1 st Monday in September	Labor Day
4 th Thursday in November	Thanksgiving
4 th Friday in November	Friday after Thanksgiving
December 25	Christmas Day

Please note that between December 24th and January 2nd, Indiana University will be open Monday through Friday for essential operations **ONLY** and will re-open for normal operations on January 2nd. If at all possible, biological specimens for submission to Indiana University should **NOT** be collected and shipped to Indiana University after the second week of December. Should it be necessary to ship blood samples for DNA extraction to Indiana University during this period, please contact the Indiana University staff before December 20th by e-mailing alzstudy@iu.edu, so that they can arrange to have staff available to process incoming samples.

Please see: https://ncrad.org/holiday_closures.html for additional information.

- Please note that courier services may observe a different set of holidays.
- Please be sure to verify shipping dates with your courier prior to any holiday.
- **Weekend/holiday delivery must be arranged in advance with NCRAD staff.**

4.0 Globally Unique Identifier (GUID)

The GUID is a subject ID that allows researchers to share data specific to a study participant, without exposing personally identifiable information. A GUID is made up of random alpha-numeric characters and does not include any PHI in the identifier. By using GUIDs in your research data, the system can associate a single research participant's genetic, imaging, and clinical assessment data even if the data was collected at different locations or throughout different studies.

To create a GUID follow these steps:

1. Create an account: <https://bricsguid.nia.nih.gov/portal/jsp/login.jsp>
2. Once you have an account, go to the GUID Tool – Create GUID
3. To open the ‘Launch GUID Tool’ you will need to have Java installed on your device
4. In order to generate a GUID, the following PHI is required ([Appendix D](#)):
 - Complete legal given (first) name of subject at birth
 - If the subject has a middle name
 - Complete legal family (last) name of subject at birth
 - Day of birth
 - Month of birth
 - Year of birth
 - Name of city/municipality in which subject was born
 - Country of birth

5.0 NCRAD Laboratory Information

5.1 Site Required Equipment

The following materials and equipment are necessary for the processing of specimens at the collection site and are to be **supplied by the local site**:

- Personal Protective Equipment: lab coat, nitrile/latex gloves, safety glasses
- Tourniquet
- Alcohol Prep Pad
- Gauze Pad
- Bandage
- Butterfly needles and hub
- Microcentrifuge tube rack
- Sharps bin and lid
- Wet Ice Bucket
- Wet ice
- Dry ice

In order to process samples consistently across all projects and ensure the highest quality samples possible, project sites must have access to the following equipment:

- Centrifuge capable of $\geq 2000 \times g$ with refrigeration to 4°C
- -80°C Freezer

In order to ship specimens, you must provide:

- Dry ice pellets

5.2 Biospecimens Sent to NCRAD

Biospecimens collected include whole blood and CSF.

*Sites are required to contact the study partner the evening prior to pharmacokinetics visits to remind study partners not to administer the study drug at home.

- Pre-dose: prior to morning drug intake
- Post-dose: between 2 and 6 hours after the morning drug intake

Please refer to the below schedule for the biospecimen collection schedule:

	PRE-DOSE		POST-DOSE		Cerebrospinal Fluid (CSF)
	Serum (Red Top)	Plasma/Buffy Coat (EDTA, Lavender Top)	Serum (Red Top)	Plasma/Buffy Coat (EDTA, Lavender Top)	
Screening	X				X
Baseline	X	X			
Week 4	X	X	X	X	
Week 8	X	X	X	X	
Week 16 (Cohorts A & B)	X	X	X	X	
Week 24	X	X	X	X	X (2A participants only)*
Week 48	X	X	X	X	
Week 72	X	X	X	X	X*
Early Termination	X	X	X	X	X**

*CSF samples should be collected at the same time of day, either morning (between 8 and 10AM) or afternoon (between 1 and 3PM) and within 2-6 hours after the morning dose.

**ET does not require CSF draw/LP although it is encouraged.

Whole blood is collected in two different collection tubes: lavender top EDTA tubes and plain red-top serum tubes. For phase 2a of the study, at Baseline through Week 72 visits, the lavender top EDTA tube is processed locally into plasma and buffy coat fractions, aliquoted, frozen at the study site, and then shipped to NCRAD. At Screening through Week 72 the plain red-top serum tube is processed locally into serum fractions, aliquoted, frozen at the study site, and then

shipped to NCRAD. On days of CSF sampling, PK samples should be drawn at the time of lumbar puncture.

CSF will be collected and aliquoted locally, frozen at the study site, and then shipped to NCRAD. The lumbar puncture should be done **after** the MRI scan **if** performed on the same day. If the lumbar puncture is performed on a separate day from the MRI and occurs **before** the MRI, then there must be at least a 3-day window between the lumbar puncture and the MRI. If the lumbar puncture is performed on a separate day from the MRI and occurs **after** the MRI, then there is no window (waiting period) between the MRI and lumbar puncture.

Consent forms must specify that any biological samples and de-identified clinical data may be shared with academic and/or industry collaborators through NCRAD. A copy of the consent form for each subject should be kept on file by the site investigator.

The **Biological Sample and Shipment Notification Form** (see [Appendix B](#)) is completed for every visit. Be sure to complete all required fields on the **Biological Sample and Shipment Notification Form**. In addition, samples not collected need to be recorded in the notes section on the Bio Sample form. Submit a copy to NCRAD with a reason provided for the omission.

Frozen samples are to be submitted according to the shipping methods outlined in [Section 10](#). Guidelines for the processing, storage location, and timing of sample collection are listed in the following tables.

5.3 Biospecimen Collection Charts

5.3.1 Biospecimen Collection for Screening Visit

Sample Type	Tube Type	Tubes Supplied in Kit	Aliquot Volume	Tubes to NCRAD	Ship
Whole blood for serum banking	Serum (Red-Top) Tube (10 mL)	1	N/A	N/A	N/A
	Serum: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	SERUM: 0.5 ml serum aliquots per 2.0 ml cryovial	Up to 11	Frozen
CSF Collection	Sterile Containers (15 ml CSF)	(1) 2.0 ml Sarstedt tube	CSF: 2.0 ml CSF per 2.0ml Sarstedt tube	1	Frozen
		26 cryovial tubes (25 orange cap, 1 blue cap)	0.5 ml CSF aliquots per 2.0 ml orange cryovial; residual volume placed in 2.0 ml cryovial with blue cap	Up to 26	
		1 yellow cap cryovial tube	1-2 ml for local lab placed in 2.0 ml cryovial with yellow cap.	0 – do not return to NCRAD	N/A

5.3.2 Biospecimen Collection for Baseline Pre-dose, Weeks 4/8/16/48 Pre-dose and Post-dose Visits

Sample Type	Tube Type	Tubes Supplied in Kit	Aliquot Volume	Tubes to NCRAD	Ship
Whole blood for serum banking	Serum (Red-Top) Tube (10 mL)	1	N/A	N/A	N/A
	Serum: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	SERUM: 0.5 ml serum aliquots per 2.0 ml cryovial	Up to 11	Frozen
Whole blood for isolation of plasma and buffy coat	EDTA (Lavender-Top) Blood Collection Tube (10 ml)	1	N/A	N/A	N/A
	Plasma: 2.0 ml cryovials with lavender cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	PLASMA: 0.5 ml plasma aliquots per 2.0 ml cryovial	Up to 11	Frozen
	Buffy Coat: 2.0 ml cryovial	1	BUFFY COAT: 0.75 ml buffy coat aliquot	1	Frozen

5.3.3 Biospecimen Collection for Weeks 24 and 72 Pre-dose and Post-dose Visit

Sample Type	Collection Tube	Tubes Supplied in Kit	Processing/ Aliquoting	Tubes to NCRAD	Ship
Whole blood for serum banking	Serum (Red-Top) Tube (10 mL)	1	N/A	N/A	N/A
	Serum: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	SERUM: 0.5 ml serum aliquots per 2.0 ml cryovial	Up to 11	Frozen
Whole blood for isolation of plasma and buffy coat	EDTA (Lavender-Top) Blood Collection Tube (10 ml)	1	N/A	N/A	N/A
	Plasma: 2.0 ml cryovials with lavender cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	PLASMA: 0.5 ml plasma aliquots per 2.0 ml cryovial	Up to 11	Frozen
	Buffy Coat: 2.0 ml cryovial	1	BUFFY COAT: 0.75 ml buffy coat aliquot	1	Frozen
CSF Collection <i>(For CSF, we are NOT collecting both pre-dose and post-dose. Only one collection per visit.)</i>	Sterile Containers (15 ml CSF)	(1) 2.0ml Sarstedt tube	CSF: 2.0 ml CSF per 2.0ml Sarstedt tube	1	Frozen
		26 cyrovial tubes (25 orange cap, 1 blue cap)	0.5 ml CSF aliquots per 2.0 ml orange cryovial; residual volume placed in 2.0 ml cryovial with blue cap	Up to 26	Frozen
		1 yellow cap cryovial tube	1-2 ml for local lab placed in 2.0 ml cryovial with yellow cap.	0 – do not return to NCRAD	N/A

5.3.4 Biospecimen Collection for Early Termination Pre-dose and Post-dose Visit

Sample Type	Collection Tube	Tubes Supplied in Kit	Processing/ Aliquoting	Tubes to NCRAD	Ship
Whole blood for serum banking	Serum (Red-Top) Tube (10 mL)	1	N/A	N/A	N/A
	Serum: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	SERUM: 0.5 ml serum aliquots per 2.0 ml cryovial	Up to 11	Frozen
Whole blood for isolation of plasma and buffy coat	EDTA (Lavender-Top) Blood Collection Tube (10 ml)	1	N/A	N/A	N/A
	Plasma: 2.0 ml cryovials with lavender cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	PLASMA: 0.5 ml plasma aliquots per 2.0 ml cryovial	Up to 11	Frozen
	Buffy Coat: 2.0 ml cryovial	1	BUFFY COAT: 0.75 ml buffy coat aliquot	1	Frozen
CSF Collection (For CSF, we are NOT collecting both pre-dose and post-dose. Only one collection per visit.)	Sterile Containers (15 ml CSF)	(1) 2.0ml Sarstedt tube	CSF: 2.0 ml CSF per 2.0ml Sarstedt tube	1	Frozen
		26 cyrovial tubes (25 orange cap, 1 blue cap)	0.5 ml CSF aliquots per 2.0 ml orange cryovial; residual volume placed in 2.0 ml cryovial with blue cap	Up to 26	Frozen
		1 yellow cap cryovial tube	1-2 ml for local lab placed in 2.0 ml cryovial with yellow cap.	0 – do not send to NCRAD	N/A

6.0 Specimen Collection Kits, Shipping Kits, and Supplies

NCRAD will provide:

- 1) Blood sample collection kits for research specimens to be stored at NCRAD, the Blood Supplemental Supply Kit, and the Frozen Shipment Kit;
- 2) CSF collection kits including Lumbar Puncture (LP) trays, and the CSF Supplemental Supply Kit;
- 3) Clinical lab supplies (with the exception of dry ice pellets and equipment supplies listed in [Section 5.1](#)). These materials include blood tubes, pipettes, LP trays (when applicable), boxes for serum/plasma/buffy coat/CSF aliquots, as well as partially completed shipping labels to send materials to NCRAD. Kit Number Labels, Site and ADCS ID Labels, Collection and Aliquot Tube Labels will all be provided by NCRAD.

Details regarding the blood and CSF kits are found in this Manual of Procedures. Collection and Aliquot Tube Labels will be pre-printed with study information specific to the type of sample being drawn. Ensure that all tubes are properly labeled during processing and at the time of shipment according to [Section 7.1](#).

6.1 Specimen Collection Kit Contents

Collection kits contain the following (for each subject) and provide the necessary supplies to collect samples from a given subject. Do not replace or supplement any of the tubes or kit components provided with your own supplies unless you have received approval from the NCRAD Study team to do so. *Please store all kits at room temperature until use.*

VIVA-MIND Blood Collection Kit Screening

Quantity	Blood Collection Kit Components
1	Serum Red Top Blood Collection Tube (10 mL)
10	Cryovial tube (2.0 mL) with red cap
1	Cryovial tube (2.0 mL) with blue cap
3	Disposable graduated transfer pipette
12	Pre-printed Aliquot Tube Label
4	Pre-printed Kit Number Label
3	Labels for Handwritten Site and ADCS ID
1	Microcentrifuge box (48-slot)
1	Resealable bag

**VIVA-MIND Blood Collection Kit- Baseline Pre-dose, Weeks 4/8/16/24/48/72, ET Pre-dose,
Post-dose**

Quantity	Blood Collection Kit Components
1	Serum Red Top Blood Collection Tube (10 mL)
1	EDTA Lavender Top Blood Collection Tube (10 mL)
10	Cryovial tube (2.0 mL) with lavender cap
1	Cryovial tube (2.0 mL) with gray cap
10	Cryovial tube (2.0 mL) with red cap
2	Cryovial tube (2.0 mL) with blue cap
3	Disposable graduated transfer pipette
23	Pre-printed Aliquot Tube Label
4	Pre-printed Kit Number Label
3	Labels for Handwritten Site and ADCS ID
1	Microcentrifuge box (48-slot)
1	Resealable bag

VIVA-MIND Blood Supplemental Supply Kit

Quantity	Blood-Based Supplemental Supply Kit Components
5	Serum (Red-Top) Blood Collection Tube (10 ml)
5	EDTA (Lavender-Top) Blood Collection Tube (10 ml)
10	Disposable graduated transfer pipette
10	Labels for handwritten Site and ADCS ID
5	Microcentrifuge box (48-slot)
10	Plastic Biohazard bag with absorbent sheet (small)
3	Warning label packet with dry ice sticker
2	Fine point permanent marker

VIVA-MIND Frozen Shipping Kit (Batch)

Quantity	Frozen Shipping Kit Components
5	Plastic Biohazard bag with absorbent sheet (small)
1	Shipping box/Styrofoam container
1	Warning label packet with dry ice sticker

VIVA-MIND LP Kits

***Sites must specify 22 or 24 gauge kit when ordering from NCRAD.**

Quantity	LP Kit Components
1	Sprotte needle, 22 or 24 gauge X 3.5" (90mm)
1	Introducer needle, 1 mm x 30 mm
1	Hypodermic needle, 22 gauge x 1.5"
1	Plastic syringe, (3 ml, luer lock) with 25G x 5/8" needle attached
4	Polypropylene syringe (5 ml, luer lock)

1	Needle stick pad
1	Adhesive bandage
1	Drape, fenestrated, 2 tabs, paper, 18" x 26"
2	Towel, 13.5" x 18"
6	Gauze pad, 2" x 2"
3	Sponge stick applicator
2	Lidocaine 1%, 5 ml
1	Povidone-Iodine Topical Solution, 0.75 oz

VIVA-MIND CSF Kits Screening

Quantity	CSF Kit Components
2	15 ml conical polypropylene tube-individually wrapped
1	50ml conical polypropylene tube
1	Bubble wrap bag
1	Large Biohazard bag with absorbent sheet
1	2.0ml Sarstedt tube
25	Cryovial tube (2.0 ml) with orange cap
1	Cryovial tube (2.0 ml) with blue cap
1	Cryovial tube (2.0 ml) with yellow cap
3	Pre-printed Kit Number labels
1	Site/ADCS ID Labels
27	Pre-printed Aliquot Tube Label
2	Disposable pipet
1	Shipping box/Styrofoam container
1	Warning label packet with dry ice sticker
1	Microcentrifuge box (48-slot)

VIVA-MIND CSF Kits W24, 72, ET

Quantity	CSF Kit Components
2	15 ml conical polypropylene tube-individually wrapped
1	50ml conical polypropylene tube
1	Bubble wrap bag
1	Large Biohazard bag with absorbent sheet
1	2.0ml Sarstedt tube
25	Cryovial tube (2.0 ml) with orange cap
1	Cryovial tube (2.0 ml) with blue cap
1	Cryovial tube (2.0 ml) with yellow cap
3	Pre-printed Kit Number labels
1	Site/ADCS ID Labels
27	Pre-printed Aliquot Tube Label

2	Disposable pipet
1	Microcentrifuge box (48-slot)

Supplemental CSF Kits

Quantity	CSF Supplemental Supply Kit Components
10	50ml conical polypropylene tube
10	15ml conical polypropylene tube-individually wrapped
5	2.0ml Sarstedt tube
5	3 ½" × 22 or 24G Sprotte needle with Introducer (90mm)

Individual Supplies

Quantities	Items Available upon request within the NCRAD kit module.
By Request	Microcentrifuge box (48-slot)
By Request	Cryovial tube (2.0 ml) with lavender cap
By Request	Cryovial tube (2.0 ml) with red cap
By Request	Cryovial tube (2.0 ml) with orange cap
By Request	Cryovial tube (2.0 ml) with yellow cap
By Request	Cryovial tube (2.0 ml) with blue cap
By Request	Cryovial tube (2.0 ml) with gray cap
By Request	2.0ml Sarstedt tube
By Request	15ml conical polypropylene tube-individually wrapped
By Request	Bubble wrap bag
By Request	50ml conical polypropylene tube
By Request	Plastic biohazard bag with absorbent sheet (large)
By Request	Plastic biohazard bag with absorbent sheet (small)
By Request	Disposable graduated transfer pipette
By Request	Plain Red Top Serum (Red-Top) Blood Collection Tube (10 ml)
By Request	EDTA (Lavender-Top) Blood Collection Tube (10 ml)
By Request	UN3373 label
By Request	Biohazard label
By Request	Dry ice shipping label
By Request	Fine Point Markers
By Request	Site and ADCS ID Labels

6.2 Kit Supply to Study Sites

Each individual site will be responsible for ordering and maintaining a steady supply of kits from NCRAD. We advise sites to keep a supply of each kit type available. Be sure to check your supplies and order additional materials before you run out or supplies expire so you are prepared for study visits. Please go to

kits.iu.edu/vivamind to request additional kits and follow the prompts to request the desired supplies. Options include ordering a specific number of kits; we are also including the option of simply ordering the desired amount of extra supplies.

Please allow **TWO weeks** for kit orders to be processed and delivered.

7.0 Blood Collection and Processing Procedures

Important Note

In order to ensure the highest quality samples are collected, processed, and stored, it is essential to follow the specific collection, processing, and shipment procedures detailed in the following pages. Please read the following instructions first before collecting any specimens. Have all your supplies and equipment out and prepared prior to drawing blood. Please note that the centrifuge may take 30 minutes to cool, so please plan accordingly. Draw blood in the following order:

1. Plain Red Top Serum Blood Collection Tube (10 ml)
2. EDTA (Lavender-Top) Blood Collection Tube (10 ml) for Buffy Coat and Plasma


7.1 Labeling Samples

Label Type Summary


1. Kit Number Label
2. Site and ADCS ID Label
3. Aliquot Tube Label
4. Sarstedt CSF Tube Label



The **Kit Number Labels** do not indicate a specimen type but are affixed on the Biological Sample and Shipment Notification Forms, each collection tube label in the kit, and on the lid of the cryobox holding the samples for the kit.



The **Site and ADCS Labels** are placed on all collection tubes, both blood and CSF.




The **Aliquot Tube Label** is placed on each blood derivative and CSF aliquot.



The **Sarstedt CSF Tube Label** is placed on the 2ml Sarstedt tube derivative.

****Important Note****


Each blood collection tube will contain two labels: the Kit Number Label and the Site and ADCS ID Label. The CSF Sarstedt Tube will only contain one label due to limited space. Be sure to place labels in the same configuration consistently among tubes, with the barcoded label near the top of the tube and the handwritten Site and ADCS label.



In order to ensure the label adheres properly and remains on the tube, please follow these instructions:

- Place cryovial labels on **ALL** cryovials **BEFORE** sample collection, processing, or freezing. This should help to ensure the label properly adheres to the tube before exposure to moisture or different temperatures.

- Using a fine point permanent marker, fill-in and place the ADCS ID label on the collection tubes **BEFORE** sample collection or processing. This label is in addition to the kit number label. **DO NOT** place ADCS ID labels on any cryovials.
- Place cryovial label **horizontally** on the tube (wrapped around sideways if the tube is upright) and just below the cap. **DO NOT** cover the barcode on the cryovial with the cryovial label (see following diagram).
- Take a moment to ensure the label is **completely adhered** to each tube. It may be helpful to roll the tube between your fingers after applying the label.



If there are any unused cryovials, please do not send the empty cryovials to NCRAD. These unused cryovials (ensure labels are removed) can be saved as part of a supplemental supply at your site or the cryovials can be disposed of per your site's requirements.

7.2 Video List

The following training videos are available to assist you with the specimen processing, aliquoting, and shipping processes. The videos are available at: <https://ncrad.org/resource/viva-mind.html>

- VIVA-MIND MOP Training
- Plasma and Buffy Coat Processing and Aliquoting
- Serum Processing and Aliquoting
- Frozen Shipping

7.3 Filling Aliquot Tubes (Plasma, Serum, and CSF)


In order to ensure that NCRAD receives a sufficient amount of sample for processing and storage, and to avoid cracking of the tubes prior to shipment, each cryovial should be filled to the assigned volume with the respective biological

material after processing is completed (refer to detailed processing instructions for average yield per sample).

Over-filled tubes may burst once placed in the freezer, resulting in a loss of that sample.

Aliquot the remaining biologic material as the residual volume and ship to NCRAD. Essentially, all material should be shipped to NCRAD, ensuring maximum amount in as many cryovials as will allow after processing the sample. For example, if 3.3 ml of sample is obtained, you should fill 6 cryovial tubes each with 0.5 ml, and one additional cryovial tube with the remaining 0.3 ml for plasma, serum, and CSF.

NOTE: 2.0ml of CSF must be drawn directly into the 2.0ml Sarstedt tube



Please note: It is critical for the integrity of the samples that study staff note if an aliquot tube contains a residual volume (anything under 0.5 ml for plasma, serum, and CSF). Please record the specimen number and volume of the residual aliquot on the Biological Sample and Notification Form.

To assist in the preparation and aliquoting of samples, colored caps and cap stickers are used for the cryovial tubes. The chart below summarizes the association between cap color and type of cryovial.

Cap Color	Sample Type
Lavender Cap	Plasma
Red Cap	Serum
Blue Cap	Residual
Gray Cap	Buffy Coat
Orange Cap	CSF
Yellow Cap	CSF for Local Lab
Blue Cap – Individually wrapped	Sarstedt tube



7.4 Serum (Red-top) Tube (10 mL) for Serum

Whole Blood Collection for Isolation of Serum: Serum (Red-Top) Tube (10 ml) (for processing of serum aliquots). One Red-Top tube is collected at Screening, Baseline, Weeks 4/8/16/24/48/72 and/or Early Termination, when applicable.

1. Set centrifuge 4°C to pre-chill before use.
2. Place completed Site and ADCS ID Label and Collection “**SERUM**” Tube Labels on the Plain Red-Top Serum Blood Collection Tube. Place pre-printed Aliquot “**SERUM**” Tube Labels on the (10) 2.0 ml cryovial tubes with red caps and (1) 2.0 ml cryovial with blue cap (if necessary, for residual).
3. Using a blood collection set and a holder, collect blood into **Plain Red-Top Serum Blood Collection Tubes (10 ml)** using your institution's recommended procedure for standard venipuncture technique.


The following techniques shall be used to prevent possible backflow:

- a. Place donor's arm in a downward position.
 - b. Hold tube in a vertical position, below the donor's arm during blood collection.
 - c. Release tourniquet as soon as blood starts to flow into tube.
 - d. Make sure tube additives do not touch the stopper or the end of the needle during venipuncture.
4. Allow at least 10 seconds for a complete blood draw to take place in each tube. **Ensure that the blood has stopped flowing into each tube before removing the tube from the holder.** The tube with its vacuum is designed to draw 5 ml of blood into the tube.
 - a. If complications arise during the blood draw, please note the difficulties on the ‘Biological Sample and Shipment Notification Form’. Do not attempt to draw an additional Serum tube at this time. Process blood obtained in existing Serum tube.
 5. **CRITICAL STEP:** Immediately after blood collection, **gently** invert/mix (180 degree turns) each tube 5 times.
 6. **CRITICAL STEP:** Allow blood to clot at room temperature by placing it upright in a vertical position in a tube rack for 30 minutes. If sample is not clotted allow it to set up to 60 minutes to clot. Serum samples need to be spun, aliquoted, and placed in the freezer within 2 hours from the time of collection.

7. After 30 minutes of clotting, centrifuge the collection tube for 10 minutes at 2000 x g at 4°C. **It is critical that the tube be centrifuged at the appropriate speed to ensure proper serum separation** (see worksheet in [Appendix A](#) to calculate RPM)
 - a. Equivalent rpm for spin at 2000 x g
 - b. While centrifuging, remember to record all times, temperatures and spin rates on the Biological Sample and Shipment Notification Form [Appendix B](#).
 - c. Serum samples need to be spun, aliquoted, and placed in the freezer within 2 hours from the time of collection.
 - d. Record time aliquoted on the Biological Sample Shipment and Notification Form. Date and time of study drug intake/dose, the time of last meal, and the time of the blood draw (pre-dose (prior to morning drug intake) and post-dose (between 2 and 6 hours after intake of the morning dose)) should be recorded on the Biological Sample Shipment and Notification Form for the day of the study visit.

8. Remove the serum by tilting the tube and placing the pipette tip along the lower side of the wall. Using a disposable pipette, transfer serum into the pre-labeled cryovials with the red caps. Aliquot 0.5 ml per cryovial (total vials = up to ten with 0.5 ml each or nine with 0.5 mL and one residual with <0.5 ml). Be sure to only place **serum** in cryovials labeled with the “SERUM” label and red caps. If there is extra serum left, use 1 extra blue-cap cryovial provided for another <0.5 ml aliquot of serum and label as appropriate. **If a residual aliquot is created, document the sample number and volume on the Biological Sample and Shipment Notification Form.**

9. Place the labeled cryovials in a cryobox and place on dry ice. Transfer to **-80°C Freezer when possible**. Store all samples at **-80°C until shipped** to NCRAD on dry ice pellets. Record time aliquots placed in freezer and storage temperature of freezer on Biological Sample and Shipment Notification Form.




Serum Aliquots: up to 10 (0.5ml each) red cap cryovials possible or nine red capped cryovials (0.5ml) and one blue capped residual (<0.5ml)

Serum Preparation (10ml Red Top Tube)




Step One




- Store tubes at room temperature.
- Label tubes with pre-printed labels prior to blood draw.

Step Two



- Collect blood in Serum Tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three




- Immediately after blood draw, invert tubes 5 times to mix samples.

Step Four



- Allow blood to clot for 30 minutes.
- Centrifuge samples at 2000 x g for 10 minutes at 4°C.
- Serum samples need to be spun, aliquoted, and in the freezer within 2 hours from the time of collection.

Step Five



- Label cryovial tubes with preprinted labels.
- Aliquot 0.5 ml into each cryovial tube.
- If residual aliquot is created, use blue cap to indicate volume difference and document Specimen Number on Biological Sample and Shipment Notification Form.
- Store serum aliquots at -80°C until shipment.

7.5 EDTA (Lavender-Top) Blood Collection Tube (10 mL) for Plasma and Buffy Coat

Whole Blood Collection for Isolation of Plasma and Buffy Coat: EDTA (Lavender-Top) Blood Collection Tube (10 ml) (for processing of plasma aliquots and buffy coat aliquot).


1. Set centrifuge to 4°C to pre-chill before use.
2. Place completed Site and ADCS ID Label and pre-printed “**PLASMA**” Collection and Tube Label on the lavender-top EDTA tube. Place pre-printed “**PLASMA**” Aliquot Tube Labels on the ten (10) 2.0 ml cryovial tubes with lavender caps. Place pre-printed “**BUFFY COAT**” Collection and Aliquot Tube Label on the (1) 2 ml cryovial with a gray cap.
3. Please ensure that aliquots are kept in numerical order (by specimen number) throughout the aliquoting and shipping process.
4. Using a blood collection set and a holder, collect blood into the **EDTA (Lavender-Top) Blood Collection Tube (10 ml)** using your institution's recommended procedure for standard venipuncture technique.

The following techniques shall be used to prevent possible backflow:

- a. Place donor's arm in a downward position.
 - b. Hold tube in a vertical position, below the donor's arm during blood collection.
 - c. Release tourniquet as soon as blood starts to flow into tube.
 - d. Make sure tube additives do not touch stopper or end of the needle during venipuncture.
5. Allow at least 10 seconds for a complete blood draw to take place in each tube. **Ensure that the blood has stopped flowing into the tube before removing the tube from the holder.** The tube with its vacuum is designed to draw 10 ml of blood into the tube.
 - a. If complications arise during the blood draw, please note the difficulties on the 'Biological Sample and Shipment Notification Form'. Do not attempt to draw an additional EDTA tube at this time. Process blood obtained in existing EDTA tube.
 6. **CRITICAL STEP:** Immediately after blood collection, **gently** invert/mix (180 degree turns) the EDTA tube 8-10 times.

7. **CRITICAL STEP:** Immediately after inverting the EDTA tube, place it on wet ice until centrifugation begins.
 - a. Preferably within 30 minutes of blood collection, centrifuge balanced tubes for 10 minutes at 2000 RCF (x g) at 4°C. **It is critical that the tubes be centrifuged at the appropriate speed and temperature to ensure proper plasma separation** (see worksheet in [Appendix A](#) to calculate RPM).
 - b. Equivalent rpm for spin at 2000 x g
 - c. While centrifuging, remember to record all times, temperatures and spin rates on the Biological Sample and Shipment Notification Form.
 - d. **Plasma samples need to be spun, aliquoted, and placed in the freezer within 2 hours from the time of collection.**
 - e. Record time aliquoted on the Biological Sample and Shipment Notification Form. Date and time of study drug intake/dose, the time of last meal, and the time of the blood draw (pre-dose (prior to morning drug intake) and post-dose (between 2 and 6 hours after intake of the morning dose)) should be recorded on the Biological Sample Shipment and Notification Form for the day of the study visit.

8. Remove the plasma, being careful not to agitate the packed red blood cells at the bottom of the collection tube. Tilt the tube and placing the disposable pipette tip along the lower side of the wall without touching the pellet (buffy coat) so that plasma is not contaminated (see below). Transfer plasma into the pre-labeled cryovials. Aliquot 0.5 ml per cryovial (up to 10 vials with 0.5 ml each). Be sure to only place **plasma** in cryovials labeled with “PLASMA” labels. Take caution not to disturb the red blood cells at the bottom of the tube. If there is extra plasma left, use 1 extra cryovial with blue cap provided for another <0.5 ml aliquot of plasma. **If a residual aliquot (<0.5 ml) is created, document the sample number and volume on the Biological Sample and Shipment Notification Form.**




NOTE: When pipetting plasma from the plasma tube into the cryovials, be very careful to pipette the plasma top layer only, leaving the buffy coat and the red blood cell layers untouched.

9. Place the labeled cryovials in the cryobox and place upright on dry ice. Transfer to **-80°C Freezer when possible**. Store all samples upright at **-80°C until shipped** to NCRAD on dry ice pellets. Record time aliquots placed in freezer and storage temperature of freezer on Biological Sample and Shipment Notification Form.

10. After plasma has been removed from the EDTA (Lavender-Top) Blood Collection Tube (10 ml), aliquot buffy coat layer (in the top layer of cells, the buffy coat is mixed with RBCs-see figure) into labeled cryovial with gray cap using a disposable graduated micropipette. All of the buffy coat will be placed into one cryovial. The buffy coat aliquot is expected to have a reddish color from the RBCs. Be sure to place buffy coat into cryovial with the gray cap and "BUFFY COAT" label.



Plasma Aliquots: up to 10 lavender capped cryovials (0.5ml each) possible or nine lavender capped (0.5ml) and one blue capped residual (<0.5ml)
 Buffy Coat Aliquots: One gray capped cryovial



11. Dispose of collection tube with red blood cell pellet according to your site's guidelines for disposing of biomedical waste.
12. Place the labeled cryovial in a cryobox and place on dry ice. Transfer to **-80°C Freezer when possible**. Store all samples upright at **-80°C until shipped** to NCRAD on dry ice pellets.

Plasma and Buffy Coat Preparation (10ml Lavender-Top Tube)



Step One




- Store tubes at room temperature.
- Label tubes with preprinted labels prior to blood draw.

Step Two



- Collect blood in Plasma Tube allowing blood to flow for 10 seconds and ensuring blood flow has stopped.

Step Three




- Immediately after blood draw, invert tubes 8-10 times to mix samples.

Step Four




- Place thoroughly mixed tube on wet ice until centrifugation begins.

Step Five



- Preferably within 30 minutes of blood draw, centrifuge samples at 2000 x g at 4°C for 10 minutes.
- Samples need to be spun, aliquoted, and in the freezer within 2 hours from the time of collection.



Step Six

- Adhere preprinted labels to the lavender cap cryovials.
- Aliquot 0.5 ml into each cryovial tube.
- If a residual aliquot is created, document specimen number and volume on Sample Notification Form.
- Store plasma aliquots at -80°C until shipment.

Step Seven

- Adhere preprinted "BUFFY COAT" labels to the gray cap cryovial.
- Using a clean pipette tip, collect the buffy coat (may have residual plasma and some RBCs included).
- Transfer the buffy coat from each EDTA tube into its own cryovial tube.
- Store buffy coat aliquot at -80°C until shipment to NCRAD.

8.0 Cerebrospinal Fluid Collection and Processing Procedures

*****Important Note*****

Fasting prior to the lumbar puncture is not required. However, efforts should be made to ensure each LP done on a participant is performed at the same time of day, either morning (between 8 and 10AM) or afternoon (between 1 and 3PM).

There are general guidelines to follow in regard to CSF Collection.

- Begin by confirming participant consented to lumbar puncture (LP) before scheduling the procedure and again prior to performing procedure.
- Do NOT use any extension tubing due to the tendency of manufactured plastic tubing to bind beta amyloid peptides and other important AD biomarkers.
- If LP was attempted but unsuccessful in obtaining CSF, a second attempt under fluoroscopy (if deemed appropriate by site clinician) is allowed.
- LP under fluoroscopy is permitted, if needed. Site personnel should advise the subject that use of fluoroscopy (x-rays) involves exposure to radiation.
- Ensure you have at least two “Lumbar Puncture Tray Kits” and sufficient “CSF Supplemental Supply Kit” provisions on hand prior to scheduling an LP visit. Also ensure adequate site-provided supplies (see above), including pelleted dry ice. Check expiration dates on all supplies, especially lidocaine.

8.1 Scheduling the LP

LPs can be performed in the morning (between 8 and 10AM) or in the afternoon (between 1 and 3PM). CSF amyloid levels can vary depending upon the time of day the sample is collected. It is important for the time of day of collection to remain consistent across study visits.


The LP should be rescheduled if the participant does not feel well or is febrile.

8.2 Performing the LP

Leave everything wrapped until the person performing the lumbar puncture is seated.

Feel the outside of the lumbar puncture kit (still wrapped up) to determine which end contains the spongy swabs. Turn this end toward the person performing the lumbar puncture and begin un-wrapping the kit.

Lumbar Puncture Tray Kit Images



Exterior of LP Tray provided by NCRAD which contains the 22 gauge Sprotte Needle with Introducer



Interior of LP Tray Provided by NCRAD

Close up of Sprotte Spinal Needle (22 gauge x 3 ½ in.) with Introducer (24 gauge is equivalent but with lavender top needle)

TOUCH ONLY THE OUTSIDE OF THE PAPER WRAPPER

When you grab an edge to unfold it, touch only the folded under portions of the outside of the wrapper. Also, don't let the outside of the wrapper touch any part of the inside.

- If you touch any part of the paper wrapper, or if any non-sterile object outside of the wrapper touches any part of the inside of the wrapper, throw the kit away and start over.
- If you are in any doubt as to whether the inside of the wrapper has been touched, throw the kit away and start over.

Unwrapping the Sterile 15 ml Conical Tubes

Note that the 15ml conical tubes and Sarstedt collection tube, into which CSF is collected and transferred come individually wrapped and are sterile inside and out. These wrappers should be peeled open by an assistant (not touching the tube) and the tube carefully dropped onto the LP tray or elsewhere in the sterile field in a manner that avoids contamination. Any additional needles or other individually wrapped sterile items can be handled the same way.

- Do not drop any packaging onto the tray or sterile field.
- Do not let the item touch the outside of the packaging on its way to the tray.

General CSF Collection Methods

Reminder: 15 ml is the required MINIMUM for CSF biomarker analysis. After the yellow capped cryovial for local lab testing has been collected, the 2ml Blue Sarstedt collection tube MUST be collected next before collecting CSF in the 15 ml conical tube. If the full 15 ml or the 2ml Blue Sarstedt collection tube is not obtained and provided to NCRAD, document the reason for under-collection and/or why the 2ml Blue Sarstedt collection tube was not collected on the comments section of the CSF Sample and Shipment Notification Form and email coordinator Dione Keys at: dlkeys@iupui.edu.

Step by Step Summary of CSF Collection Procedure

1. Print CSF Sample and Shipment Notification Form.
2. Confirm all supplies, including dry ice **pellets** and wet ice, are available.
3. Label the (25) orange cap cryovials, and (1) blue cap cryovial with provided VIVA-MIND CSF labels. Do **NOT** open and label the individually wrapped 15 ml conicals or Sarstedt tube that will be kept sterile to collect the CSF.
4. Pre-cool the centrifuge and pre-cool all labeled cryovials on wet ice. Do **NOT** pre-cool the 15 ml conicals or Sarstedt tube that will be kept sterile to collect the CSF.

5. Measure vitals (participant lying down).
6. Record the time of LP and associated information on the CSF Sample and Shipment Notification Form.
7. Collect CSF (by gravity-per protocol) and document method on CSF Sample and Shipment Notification Form) following these steps:
 - a. Collect initial 1-2 ml (if bloody, collect CSF until cleared of blood) using the 15ml conical tube. If not bloody, transfer first 1-2ml into yellow cap cryovial for local lab.
 - b. **Next, collect 2.0 mls directly into the 2.0ml Sarstedt tube. The tube has a “max fill” line at the 2.0 ml mark. No processing is allowed for this tube – at no point is the tube centrifuged or aliquoted to/from. Label the (1) 2.0ml Sarstedt tube, cap the tube and store.**
 - a. **** 2ml Sarstedt tube MUST be collected before aliquoting into the orange cap cryovials.**
 - c. Collect additional CSF (totaling 15ml including the 2.0ml-3.0ml collected initially) into the **UNLABELED-STERILE** 15ml polypropylene conicals from the “CSF Supply Kit”.
8. Process the CSF as follows:
 - a. Place samples upright on wet ice prior to processing. Preferably within 15 minutes of collection, centrifuge briefly at low speed (2000 x g, 10 min, 4°C with brake turned off) to pellet any cellular debris.
 - b. Aliquot 0.5ml into the orange-cap cryovials. If a residual aliquot is created, aliquot into blue-cap cryovial. Document specimen number and volume on CSF Sample and Shipment Notification Form.
 - c. Store CSF aliquots at -80°C and record time of freezing on CSF Sample and Shipment Notification Form.
9. **REMINDER:** If the full 15 ml or the 2ml Blue Sarstedt collection tube is not obtained and provided to NCRAD, document the reason for under-collection and/or why the 2ml Blue Sarstedt collection tube was not collected on the comments section of the CSF Sample and Shipment Notification Form and email coordinator Dione Keys at: dlkeys@iupui.edu.

CSF Preparation (15 ml total)


****2ml Sarstedt tube MUST be collected**

Step One



- Label tubes with pre-printed subject labels prior to collection.
- Pre-chill all cryovials on wet ice.

Step Two



- Collect initial 1-2ml (if bloody, collect CSF until cleared of blood) into 15 ml conical tube.
- If not bloody, transfer 1-2ml into the yellow cap cryovial and send to local lab for testing.

Step Three



- Next collect 2.0 ml into the 2.0ml Sarstedt tube.
- The tube has a "max fill" line at the 2.0ml mark.
- No processing is allowed for the Sarstedt tube.
- Label the Sarstedt tube, cap the tube and store.

Step Four




- Collect 15 ml total, including the 2.0-4.0 ml collected in steps two and three.

Step Five




- Place 15ml conical upright on wet ice until centrifugation begins.

Step Six



- Preferably within 15 minutes of collection, centrifuge samples at 4°C at 2000 x g for 10 minutes.

Step Seven



- Using a clean pipette, aliquot 0.5 ml into the orange-cap cryovials, leaving the debris pellet undisturbed.
- If a residual aliquot is created, aliquot into blue-cap cryovial. Document specimen number and volume on CSF Sample Notification Form.
- Within 2 hours of CSF collection, samples need to be spun, aliquoted and in the freezer. Store at -80°C until shipment. Record time of freezing on CSF Sample Notification Form.

9.0 Sample Redraws

*****Important Note*****

If challenges arise during the blood draw process, it is advised that the phlebotomist discontinue the draw. Attempt to process and submit any blood-based specimens that have already been collected to NCRAD.

Redraws will be scheduled for samples submitted to NCRAD.

There may be situations that arise that require a patient sample to be redrawn from certain cycles/visits. At those times, NCRAD study staff will alert site coordinators that a participant sample has failed and should be redrawn. This can happen for several reasons, including insufficient blood at the time the sample was drawn, temperature storage extremes, or even shipping errors.

1. If the biospecimens at a scheduled visit **are partially** collected:
 - a. Attempt to process and submit any samples that were able to be collected during the visit.
 - b. Document difficulties on the 'Biological Sample and Shipment Notification Form' prior to submission to NCRAD.
 - i. Indicate blood draw difficulties at the bottom of the 'Biological Sample and Shipment Notification Form' within the "Notes" section.
 - ii. Complete the 'Biological Sample and Shipment Notification Form' with tube volume approximations and number of aliquots created.
 - c. Contact a NCRAD coordinator and alert them of the challenging blood draw.

2. If the biospecimens at a scheduled visit **are not** collected:
 - a. Contact the VIVA-MIND Project Manager and a NCRAD coordinator to alert them of the challenging blood draw or circumstances as to why biospecimens were not collected.

Schedule participant for a re-draw visit as quickly as possible.

10.0 Packaging and Shipping Instructions

ALL study personnel responsible for shipping should be certified in biospecimen shipping. If not available at your University, please contact NCRAD with questions and information regarding resources.

Sample Type	Collection Tube	Tubes Supplied in Kit	Processing/ Aliquoting	Tubes to NCRAD	Ship
Whole blood for serum banking	Serum (Red-Top) Tube (10 mL)	1	N/A	N/A	N/A
	Serum: 2.0 ml cryovials with red cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	SERUM: 0.5 ml serum aliquots per 2.0 ml cryovial	Up to 11	Frozen
Whole blood for isolation of plasma and buffy coat	EDTA (Lavender-Top) Blood Collection Tube (10 ml)	1	N/A	N/A	N/A
	Plasma: 2.0 ml cryovials with lavender cap (residual volume placed in 2.0 ml cryovial with blue cap)	11	PLASMA: 0.5 ml plasma aliquots per 2.0 ml cryovial	Up to 11	Frozen
	Buffy Coat: 2.0 ml cryovial	1	BUFFY COAT: 0.75 ml buffy coat aliquot	1	Frozen
CSF	Sterile Containers (15 ml CSF)	27	CSF: 2.0 ml CSF per 2.0ml Sarstedt tube	1	Frozen
			CSF: 0.5 ml CSF aliquots per 2.0 ml orange cryovial	Up to 26	

Specimens being shipped to NCRAD should be considered as Category B UN3373 specimens and as such must be tripled packaged and compliant with IATA Packing Instructions 650. *See the Latest Edition of the IATA Regulations for complete documentation.*

Triple packaging consists of a primary receptacle(s), a secondary packaging, and a rigid outer packaging. The primary receptacles must be packed in secondary packaging in such a way that, under normal conditions of transport, they cannot break, be punctured, or leak their contents into the secondary packaging. Secondary packaging must be secured in outer packaging with suitable cushioning material. Any leakage of the contents must not compromise the integrity of the cushioning material or of the outer packaging.

10.1 Frozen Shipping Instructions

*****Important Note*****

Screening CSF samples should be shipped every week, Monday-Wednesday to NCRAD.

All blood samples, and Weeks 24, 48, 72 and Early Termination CSF samples should be batch shipped to NCRAD in increments of five (5) cryoboxes, or every three (3) months, whichever comes first.

The most important issue for shipping is to maintain the temperature of the samples. The frozen samples must never thaw; not even the outside of the tubes should be allowed to defrost. **This is best accomplished by making sure the Styrofoam container is filled completely with pelleted dry ice.**



Large Frozen Shipper:

**** 45 lbs of dry ice pellets**

AND

- Fits up to 5 48-slot cryoboxes



Small Frozen Shipper:

**** 10 lbs of dry ice pellets**

AND


- fits up to 2 x 48-slot cryoboxes
- OR**
- fits 1 x 48-slot cryobox and 2ml Sarstedt Tube
- OR**
- 1 x 2ml Sarstedt tube

*****Important Note*****

FROZEN SAMPLES MUST BE SHIPPED MONDAY-WEDNESDAY ONLY!

***** Packing and Labeling Guidelines *****


- The primary receptacle (cryovial) must be leak proof and must not contain more than 1L total.
- The secondary packaging (biohazard bag) must be leak proof and if multiple blood tubes are placed in a single secondary packaging, they must be either individually wrapped or separated to prevent direct contact with adjacent blood tubes.
- Absorbent material must be placed between the primary receptacle and the secondary packaging. The absorbent material should be of sufficient quantity in order to absorb the entire contents of the specimens being shipped. Examples of absorbent material are paper towels, absorbent pads, cotton balls, or cellulose wadding.
- A shipping manifest of specimens being shipped must be included between the secondary and outer packaging.
- The outer shipping container must display the following labels:
 - ✓ Sender's name and address
 - ✓ Recipient's name and address
 - ✓ Responsible Person
 - ✓ The words "Biological Substance, Category B"
 - ✓ UN3373
 - ✓ UPS Dry Ice label and net weight of dry ice contained



Packaging Instructions


1. Contact UPS to confirm that service is available and arrange for the package to be picked up.
2. Notify NCRAD of shipment by emailing NCRAD coordinators at alzstudy@iu.edu. Attach the following to the email:
 - a. Completed Sample Forms ([Appendix B](#) and [Appendix C](#)) to the email notification (email NCRAD coordinator prior to shipment to receive sample form).
 - b. If email is unavailable, please call NCRAD at 1-800-526-2839 and do not ship until you've contacted and notified NCRAD coordinators about the shipment in advance.
3. Place the cryovial boxes containing frozen labeled samples into a biohazard bag.

- a. **BLOOD KITS:** Place plasma, buffy coat and serum samples within one cryobox per participant.




Place kit number label on each cryobox

- b. **CSF:** Place CSF samples within one cryobox per participant.




Place kit number label on each cryobox

- c. **SCREENING VISIT ONLY:** Place 2.0ml Sarstedt tube in a 50ml conical. Place 50ml conical in bubble wrap tube sleeve. Screening visit will also have serum and CSF aliquots within one cryobox. For all other visits, CSF will be in its own cryobox.



Place kit number labels (blood and CSF) on each cryobox

- As the samples are placed in the plastic biohazard bag, do NOT remove the absorbent material found in the bag. Seal according to the instructions on the bag.



- Place approximately 2-3 inches of **dry ice pellets** in the bottom of the Styrofoam shipping container.
- Place the biohazard bags into the provided Styrofoam-lined shipping container on top of the **dry ice pellets**. Please ensure that cryovial boxes are placed so the cryovials are upright in the shipping container.
- Cover the biohazard bags containing the cryovial boxes with **dry ice pellets** to completely fill the frozen shipper.



8. After the samples have been placed into the shipping container, fill the inner Styrofoam with plenty of **dry ice pellets** to ensure the frozen state of the specimens during transit.
9. Replace the lid on the Styrofoam carton. Place the completed Blood Sample and Shipment Notification Form in the package on top of the Styrofoam lid for each patient specimen, and close and seal the outer cardboard shipping carton with packing tape.
10. Complete the UPS Dry Ice Label with the following information:
 - a. Net weight of dry ice in kg (must match amount on the airbill)
 - b. Do not cover any part of this label with other stickers, including preprinted address labels.
11. Apply all provided warning labels and the completed UPS return airbill to the outside of package, taking care not to overlap labels.
12. Hold packaged samples in -80°C freezer until time of UPS pick-up/drop-off.
13. Specimens should be sent to the following address via UPS Next Day Air. Frozen shipments should be sent Monday through Wednesday to avoid shipping delays on Thursday or Friday.

VIVA-MIND at NCRAD
Indiana University School of Medicine
351 W. 10th St. TK-217
Indianapolis, IN 46202

14. Use UPS tracking to ensure the delivery occurs as scheduled and is received by NCRAD. Please notify NCRAD by email (alzstudy@iu.edu) that a shipment has been sent and include the UPS tracking number in your email.

Shipping Instructions


*****Important Note*****

Screening CSF samples should be shipped every week, Monday-Wednesday to NCRAD.


All blood samples, and Weeks 24,48,72 and Early Termination CSF samples should be batch shipped to NCRAD in increments of five (5) cryoboxes, or every three (3) months, whichever comes first.

1. Log into the ShipExec Thin Client at kits.iu.edu/UPS.
 - a. If a new user or contact needs access, please reach out to your study contact for access.

2. Click “Shipping” at the top of the page and select “Shipping and Rating”.



3. Select your study from the “Study Group” drop down on the right side of the main screen. Choosing your study will automatically filter the address book to only addresses within this study.
4. Click on the magnifying glass icon in the “Ship From” section to search for your shipping address.



- a. Search by Company (site), Contact (name), or Address 1 (first line of your site’s street address). Click Search.
 - b. Click Select to the left of the correct contact information.
5. Verify that both the shipping information AND study reference are correct for this shipment.
 - a. If wrong study contact or study reference, click Reset in the bottom right of the screen to research for the correct information.
6. Enter Package Information
 - a. Ambient shipments

- i. Enter the total weight of your package in the “Weight” field and leave the “Dry Ice Weight” field empty.
 - b. Frozen shipments
 - i. Enter the total weight of your package in the “Weight” field.
 - ii. Enter the dry ice weight in the “Dry Ice Weight” field.
 - iii. If the “Dry Ice Weight” field is higher than the “Weight” field, you will receive an error message after clicking Ship and need to reenter these values.
 - c. Click Ship in the bottom right of the page when complete.
7. If your site does not already have a daily UPS pickup, you can schedule one here.
 - a. Click the blue Pickup Request button. Enter the earliest pickup time and latest pickup time in 24-hr format.
 - b. Give a name & phone number of someone who the UPS driver can call if having issues finding the package
 - c. Give the Floor and Room Number (if needed) to be as descriptive as possible where this package needs to be picked up from. Click Save.
8. Print the airbill that is automatically downloaded.
 - a. To reprint airbill, click History at the top left of the page.
 - i. Shipments created from the user that day will automatically populate. If shipments from a previous day need to be located, search by ship date.
 - ii. Locate the correct shipment and click on the printer icon to the left of the tracking number under “Action” to reprint the airbill.
 - iii. Click print icon on right side of the tracking number line.
9. Fold airbill and place inside plastic UPS sleeve.
10. Peel the back off of the UPS sleeve and stick the sleeve to the package.

11.0 Data Queries and Sample Reconciliation

The Laboratory worksheets must be completed on the day that samples are collected since they capture information related to the details of the sample collection and processing. These forms include information that will be used to reconcile sample collection and receipt, as well as information essential to future analyses.

Data queries or discrepancies with samples shipped and received at NCRAD may result from:

- Missing samples
- Incorrect samples collected and shipped
- Damaged or incorrectly prepared samples
- Unlabeled samples, samples labeled with incomplete information, or mislabeled samples
- Discrepant information documented on the Biological Sample and Shipment Notification Form and logged at NCRAD compared to information entered into the ADCS database.
- Samples that are frozen and stored longer than one quarter at the site
- Use of an incorrect Biological or CSF Sample and Shipment Notification Form

12.0 Appendices

Appendix A. Rate of Centrifuge Worksheet

Appendix B. Biological Sample and Shipment Notification Form

Appendix C. CSF Sample and Shipment Notification Form

Appendix A. Rate of Centrifuge Worksheet

Rate of Centrifuge Worksheet

Please complete and return this form by email to the NCRAD Project Manager if you have any questions regarding sample processing. The correct RPM will be sent back to you.

Submitter Information

Name:

Site:

Submitter e-mail:

Centrifuge Information

Please answer the following questions about your centrifuge.

Centrifuge Type

Fixed Angle Rotor:

Swing Bucket Rotor:

Radius of Rotation (mm):

Determine the centrifuge's radius of rotation (in mm) by measuring distance from the center of the centrifuge spindle to the bottom of the device when inserted into the rotor (if measuring a swing bucket rotor, measure to the middle of the bucket).

Calculating RPM from G-Force:

$$\text{RCF} = \left(\frac{\text{RPM}}{1,000} \right)^2 \times r \times 1.118 \Rightarrow \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1.118}} \times 1,000$$

RCF = Relative Centrifugal Force (G-Force)

RPM = Rotational Speed (revolutions per minute)

R= Centrifugal radius in mm = distance from the center of the turning axis to the bottom of centrifuge

Comments:

Please send this form to NCRAD Study Coordinator

alzstudy@iu.edu

Appendix D. GUID Demographics Form

Please be certain to collect the following demographic information to generate a Global Unique Identifier. This information will not be sent to NCRAD, but rather maintained at the enrolling site only:

1. Complete legal given (first) name of subject at birth: _____
2. Complete additional (middle) name or names at birth: _____
3. Complete legal family (last) name of subject at birth: _____
4. Suffix: _____
5. Date of Birth: _____
6. Name of city/municipality in which subject was born: _____
7. Country of birth: _____